
# ARROW® EZ-IO®

INTRAOSSUEOUS VASCULAR ACCESS

## EZ-IO® INTRAOSSEOUS VASCULAR ACCESS NEEDLES

### Instructions for Use

English

Teleflex®

**INDICATIONS FOR USE:**

For intraosseous access anytime in which vascular access is difficult to obtain in emergent, urgent or medically necessary cases.

**INSERTION SITES:****ADULTS**

- Proximal humerus
- Proximal tibia
- Distal tibia

**PEDIATRICS**

- Proximal humerus
- Proximal tibia
- Distal tibia
- Distal femur

**CONTRAINDICATIONS FOR USE:**

- Fracture in target bone.
- Previous, significant orthopedic procedure at the site, prosthetic limb or joint.
- IO access (or attempted IO access) in targeted bone within past 48 hours.
- Infection at the area of insertion.
- Excessive tissue (severe obesity) and/or absence of adequate anatomical landmarks.

**MRI SAFETY INFORMATION****MR Conditional**

Non-clinical testing has demonstrated that Teleflex's Arrow® EZ-IO® Needle Set is MR Conditional. A patient with these devices can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5-Tesla (1.5 T) or 3-Tesla (3 T).
- Maximum spatial field gradient of 4,000 G/cm (40 T/m).
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2.0 W/kg (Normal Operating Mode).

**RF HEATING**

Under the scan conditions defined above, Teleflex's Arrow® EZ-IO® Needle Set™ Needles are expected to produce a maximum temperature rise less than or equal to 5.1 °C after 15 minutes of continuous scanning.

**MR Artifact**

In non-clinical testing, the image artifact caused by Teleflex's Arrow® EZ-IO® Needle Set, extends approximately 6.3 cm from the device when imaged with a spin-echo or gradient-echo pulse sequence in a 3 T MRI system.

**WARNINGS AND PRECAUTIONS FOR****EZ-IO® INTRAOSSEOUS VASCULAR ACCESS SYSTEM:****CAUTIONS:**

- Use aseptic technique.
- Check skin, adipose and muscle thickness before insertion.

- Extra care should be taken during insertion and site monitoring when used in patients with bone diseases that increase the likelihood of fracture, extravasation and dislodgement.
- Do not recap Needle Sets or reconnect separated components. Use biohazard and sharps disposal precautions. Re-use of contents may cause cross-contamination, leading to patient risk and complication(s).
- Before administering vesicant, toxic, or highly-concentrated drugs, check the IO Cannula again for placement and patency.
- Use caution with chemotherapeutic agents.
- Monitor IO site/limb/infusion frequently for any signs of extravasation/infiltration, localized inflammation, changes in infusion rates or dislodgement, particularly in the first half hour after insertion, anytime the IO Cannula is manipulated or after patient transport, and during infusion of vasopressors, vesicants, and bolus or with high infusion rates and high pressure, but at least hourly during all infusions. This is especially important for all high-risk patients (elderly, pediatric, patients in shock, coagulopathies, decreased immunity, obese, etc).
- Post-IO cannula removal, a delayed complication can occur. Instruct patients and caregivers to return patient to the hospital for any problems in the limb to include a change in the limb appearance (discoloration, swelling), pain, warmth, paresthesias, fever, and prolonged discomfort.
- Complications for individuals with co-morbidities that increase risk of infection or other IO access related complications may be at a higher rate than in patients lacking co-morbidities. This risk may increase with a longer dwell/ time the device is in place.
- Do not leave the Cannula inserted for longer than 72 hours.
- Needle Sets are single use only; serious medical consequences (e.g. life-threatening infection) and reduced performance (e.g. blunted needles) may occur if compliance to this warning is not followed.
- Read all warnings, precautions, and instructions prior to use. Failure to follow these instructions and associated clinical educational materials may lead to patient or provider injury or death.
- IO infusion pain varies from mild to severe. Pain may be mitigated with a slow infusion of preservative-free and epinephrine-free lidocaine, before initial flush; and other analgesics appropriate to each patient's clinical situation.
- Potential side effects include pain, inflammation, bleeding at the insertion site, extravasation, infiltration, infection, osteomyelitis, compartment syndrome.

#### **EZ-IO® NEEDLE SETS: DESCRIPTION**

- Comprised of Cannula with Luer-lock connection, Stylet, Safety Cap.
- 15 gauge, 304 stainless steel in 15 mm, 25 mm and 45 mm lengths.
- Sterile, non-pyrogenic, in protective packaging.
- Intended for use with EZ-IO® Power Driver.

## EZ-IO® Power Driver and EZ-IO® Needle Sets: Description

EZ-IO® 45 mm Needle Set:  $\geq 40$  kg

EZ-IO® 25 mm Needle Set:  $\geq 3$  kg


EZ-IO® 15 mm Needle Set: 3-39 kg

### Insertion Instructions

For additional clinical educational resources please visit [Teleflex.com/EZIOeducation](http://Teleflex.com/EZIOeducation)

1. Clean insertion site per institutional protocol/policy.
2. Prepare supplies.
  - a. Prime EZ-Connect® Extension Set.
    - Unlock clamp.
    - Prime set and purge air.
  - b. Open EZ-Stabilizer™ Dressing package.
3. Attach EZ-IO® Needle Set to EZ-IO® Power Driver and remove Safety Cap from Cannula.  
**IMPORTANT:** Only handle EZ-IO® Needle Set by the plastic Hub.  
**IMPORTANT:** Control patient movement prior to and during procedure.
4. Push EZ-IO® Needle Set through skin until tip touches bone. 5 mm of the Cannula (at least one black line) must be visible outside the skin.

**IMPORTANT: The most accurate determinant of correct needle selection is use of depth markings. Black depth marks on each cannula function as depth measuring guides to determine soft tissue depth overlying bone (see above).**

5. Squeeze trigger and apply gentle, steady pressure.

**IMPORTANT: DO NOT USE EXCESSIVE FORCE.**

**Note:** If EZ-IO® Power Driver stalls and EZ-IO® Needle Set will not penetrate the bone, operator may be applying too much downward pressure to penetrate bone.

**Note:** In the event of an EZ-IO® Power Driver failure, disconnect the EZ-IO® Power Driver, grasp the EZ-IO® Needle Set Hub by hand and advance into the medullary space while twisting back and forth.

6. Advance EZ-IO® Needle Set and release Trigger.

**Pediatrics:** Release Trigger when sudden “give” or “pop” is felt, indicating entry into medullary space.

**Adults:** Advance EZ-IO® Needle Set approximately 1 cm after entry into medullary space; in proximal humerus for most adults Cannula should be advanced until Needle Hub is flush or against the skin (this may be more than approximately 1 cm).

7. Stabilize Needle Set Hub, disconnect EZ-IO® Power Driver, and remove Stylet.  
8. Place Stylet into NeedleVISE® for sharps containment.

**Note:** Place the NeedleVISE® on a flat stable surface. Immediately following use of a needle, use a one-handed technique holding the stylet hub, firmly insert the sharp pointed tip straight down into the opening in the NeedleVISE® until it stops. Do not hold NeedleVISE® with free hand. Dispose of opened sharp into NeedleVISE® whether or not it has been used.

9. Obtain samples for lab analysis, if needed.

**Note:** Only attach a Syringe directly to the EZ-IO® Cannula Hub when drawing blood for laboratory analysis (stabilize Cannula) or removal.

10. Place EZ-Stabilizer™ Dressing over Cannula Hub.

**Note:** Use of the EZ-Stabilizer™ Dressing is strongly recommended for all EZ-IO® Needle insertions.

11. For patients responsive to pain, consider 2% preservative-free and epinephrine-free lidocaine (intravenous lidocaine), follow institutional protocols/policy.
  - a. Local anesthetics intended for the medullary space must be administered very slowly until desired anesthetic effect is achieved.
12. Attach a primed EZ-Connect® Extension Set to the Hub, firmly secure to Cannula Hub by twisting clockwise, ensure clamp is open.

**Note:** Do NOT use any instruments to tighten connections.


**Note:** To prevent valve damage, Do NOT use needles or blunt cannula to access the swabable valve. Non-standard syringes or connectors can damage the swabable valve.

**Note:** Operator may use a sterile alcohol wipe, to swab the EZ-Connect® Extension Set valve and let it air dry.
13. Attach EZ-Stabilizer™ Dressing by pulling the tabs to expose the adhesive and adhere to skin. Secure the affected limb to minimize movement and risk of dislodgement; ambulation is discouraged. Use caution moving patients.
  - a. Proximal humerus: Secure arm in adducted position (with the patient's arm close to body), or across the abdomen using immobilizer or alternate method.
  - b. Distal Femur: Stabilize extremity and secure site with leg outstretched to ensure knee does not bend using leg board or alternate method.
  - c. Proximal and distal tibia: Minimize potential for cannula movement when necessary with use of leg board or alternate method in pediatric patients.
14. Flush the EZ-IO® Cannula with normal saline (0.9% Sodium Chloride)(5-10 mL for adults; 2-5 mL infant/child).
  - a. Prior to flush, aspirate slightly for visual confirmation of bone marrow.
  - b. Failure to appropriately flush the EZ-IO® Cannula may result in limited or no flow. Repeat flush as needed.
  - c. Once EZ-IO® Cannula has been flushed, administer fluids or medications as indicated.
15. Confirm Cannula placement with the following recommended methods:
  - Stability of Cannula in the bone.
  - Ability to aspirate after flush.
  - Adequate flow rate.
16. Document date/time of insertion and apply wristband.

**CAUTION:** Monitor insertion site frequently for extravasation.

**To remove EZ-IO® from patient:**

- Remove EZ-Connect® Extension Set.
- Lift & remove EZ-Stabilizer™ Dressing.
- Attach Luer-lock Syringe to Hub of Cannula. Maintain axial alignment and rotate clockwise while pulling straight out. Do NOT rock or bend the Cannula. Improper technique may cause cannula to break.
- Once removed, immediately place Syringe/Cannula in appropriate sharps container.
- Dress site per institutional protocol/policy.

**Note:** If the Cannula or Needle set breaks during or after placement in the patient, attempt to grasp the cannula that remains in the patient with a hemostat and remove by gently pulling while simultaneously rotating. If broken cannula is not accessible, obtain X-Ray and have physician determine if and how it should be removed as a foreign body.

Education and training materials available at [ArrowEZIO.com](http://ArrowEZIO.com)

**en** For Instructions for Use visit: [www.teleflex.com/IFU](http://www.teleflex.com/IFU)

**Teleflex®**

Customer Service: 1.866.479.8500

Manufactured for:  
Teleflex Medical  
IDA Business & Technology Park,  
Dublin Rd, Athlone,  
Co. Westmeath, Ireland

Medical  
Device

No sternal use

MR conditional

Not made with  
natural rubber latex

Sterilized using  
ethylene oxide

Do not use if  
package is  
damaged

Do not reuse

Do not  
resterilize

Keep away  
from sunlight

Keep dry

Caution

Single sterile  
barrier system

Rx only.

The System Conforms  
to the Medical Device  
Directive (93/42/EEC)

Consult  
instructions  
for use

©2021 all rights reserved. Arrow®, EZ-IO Intraosseous Vascular System®,  
EZ-Stabilizer™, and EZ-Connect® are trademarks of Teleflex Inc.

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EZ-Stabilizer™ och EZ-Connect® är varumärken som tillhör Teleflex Inc.