Sheath Adapter Product with Cath-Gard® Catheter Contamination Shield

for use with 4 Fr. - 7 Fr. Catheters

Rx only

Intended Purpose:

The Arrow® Sheath Adapter is utilized in conjunction with one-piece Percutaneous Sheath Introducer to permit venous access and catheter introduction into the central circulation.

Patient Target Group:

Intended to be used in patients with anatomy suitable for use with the introducer device.

Contraindications:

None known.

Clinical Benefits to be Expected:

Use of the Cath-Gard helps protect indwelling catheters during use from unintentional contamination.

Performance Characteristics:

Provides catheter protection from contamination.

LATEX Caution: This Product Contains Natural Rubber Latex Which May Cause Allergic Reactions.

⚠️ General Warnings and Precautions

Warnings:

1. Sterile, Single use: Do not reuse, reprocess or resterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death. Reprocessing of medical devices intended for single use only may result in degraded performance or a loss of functionality.
2. Read all package insert warnings, precautions and instructions prior to use. Failure to do so may result in severe patient injury or death.
3. Air embolism can occur if air is allowed to enter a vascular access device or vein. Do not leave open needles or uncapped, unclamped devices in venous puncture site. Use only securely tightened Luer-Lock connections with any vascular access device to guard against inadvertent disconnection. To reduce the possibility of disconnection, use bayonet lock connections only. Follow hospital protocol for sheath and adapter maintenance to guard against air embolism.
4. Sheath adapter should be removed immediately following catheter removal and replaced with Arrow obturator or obturator cap. During catheter introduction or exchange, the adapter must be occluded at all times to minimize the risk of air embolism or hemorrhage. If catheter introduction is delayed or catheter is removed, temporarily cover adapter opening with sterile-gloved finger until either catheter is inserted, sheath adapter is completely closed or an obturator or obturator cap is attached.
5. Percutaneous sheath introducers and adapters should be maintained in closely monitored environments (e.g. critical care unit, operating room, or recovery room), and not in a general nursing unit.

6. Do not apply excessive force in removing catheter or sheath. Excessive force can cause component damage or breakage. If damage is suspected or withdrawal cannot be easily accomplished, radiographic visualization should be obtained and further consultation requested.

Precautions:

1. Do not alter components during insertion, use or removal.
2. Procedure must be performed by trained personnel well versed in anatomical landmarks, safe technique and potential complications.
3. Use standard precautions and follow institutional policies for all procedures including safe disposal of devices.
4. If the package is damaged or unintentionally opened before use do not use the device. Dispose of the device.
5. Storage conditions for these devices require that they are kept dry and out of direct sunlight.
6. Sheath adapter with Cath-Gard should be used with one-piece 8.5 Fr. and 9 Fr. percutaneous sheath introducers only.
7. Perform routine insertion site care and maintenance to cap of adapter especially before inserting catheter. Make sure cap is cleaned per hospital protocol before reinserting catheter.

A Suggested Procedure: Use sterile technique.

1. Place patient in desired position.
2. Prep and drape puncture site as required.
3. Insert percutaneous sheath introducer into the patient according to the manufacturer's instructions.
4. Attach sheath adapter to hub of one-piece percutaneous sheath introducer and lock it in place via bayonet lock mechanism (refer to Figure 1). Orient the slot in sheath adapter with the locking pin on sheath hub and slide the adapter forward over sheath hub. Twist to lock.

⚠️ Warning: The connection between sheath adapter and sheath must be tightly secured and examined routinely to minimize the risk of disconnection and possible air embolism, blood loss, or exsanguination.

⚠️ Precaution: Sheath adapter should be used with one-piece 8.5 Fr. and 9 Fr. percutaneous sheath introducers only.

⚠️ Precaution: If excessive force is needed in securing adapter on sheath, do not use that adapter.


Figure 1

5. Insert tip of desired catheter through Tuohy-Borst (T-B) cap end of Cath-Gard. Advance catheter through tubing and soft hub at distal end of Cath-Gard (refer to Figure 2).

⚠️ Precaution: Sheath adapter allows passage of 4 Fr. to 7 Fr. catheters only.


Figure 2

- Feed catheter through sheath adapter and sheath into vessel and advance catheter to desired position.

Warning: Sheath adapter must be occluded at all times to minimize the risk of air embolism or hemorrhage. If catheter introduction is delayed or catheter is removed, temporarily cover adapter opening with sterile-gloved finger until either catheter is inserted, sheath adapter is completely closed or an obturator or obturator cap is attached.

- Twist sheath adapter cap clockwise to lock catheter in position. Tighten cap to the point that a slight resistance is encountered when catheter is gently tugged to ensure tight seal around it.
- Slide entire catheter contamination shield to proximal end of catheter so that soft hub of Cath-Gard is approximately 5 inches (12.7 cm) from sheath adapter (refer to Figure 3).


Figure 3

- Secure T-B cap end of Cath-Gard in place by turning the cap clockwise until slight resistance is encountered. Disengage soft hub of Cath-Gard from inner feed tube by pulling forward and advance it forward towards sheath adapter.
- Secure soft hub over adapter cap and make sure connection is secure.
- Grasp catheter through front portion of Cath-Gard catheter contamination shield and hold in place while repositioning T-B cap end as desired (refer to Figure 4). Tighten

the cap of T-B adapter also to minimize the risk of migration of Cath-Gard catheter contamination shield. Apply sterile sealing tape to T-B cap and catheter junction to establish and maintain sterile field within Cath-Gard.


Figure 4

- If catheter introduction is delayed or catheter is removed, temporarily cover adapter opening with sterile-gloved finger until either catheter is inserted or sheath adapter is completely removed.

Warning: Sheath adapter must be occluded at all times to minimize the risk of air embolism or hemorrhage.

- Remove adapter from percutaneous sheath introducer if it is no longer needed to introduce catheter and place an Arrow obturator cap or obturator via bayonet lock mechanism on the sheath hub to establish hemostasis.
- Dress sheath insertion site per hospital protocol.

Precaution: Maintain the insertion site with regular, meticulous redressing using aseptic technique.

For reference literature concerning patient assessment, clinician education, insertion technique, and potential complications associated with this procedure, consult standard textbooks, medical literature, and Arrow International LLC website: www.teleflex.com

A pdf copy of this IFU is located at www.teleflex.com/IFU

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority. The contacts of national competent authorities (Vigilance Contact Points) and further information can be found on the following European Commission website: https://ec.europa.eu/growth/sectors/medical-devices/contacts_en

en

Symbol Glossary: Symbols are in compliance with ISO 15223-1. Some symbols may not apply to this product. Refer to product labeling for symbols that apply specifically to this product.

								
Caution	Medical device	Consult instructions for use	Do not reuse	Do not resterilize	Sterilized by ethylene oxide	Single sterile barrier system	Keep away from sunlight	Keep dry
								
Do not use if package is damaged	Caution: This Product Contains Natural Rubber Latex Which May Cause Allergic Reactions.	Catalogue number	Lot number	Use by	Manufacturer	Date of manufacture	Importer	

Arrow, the Arrow logo, Cath-Gard, Teleflex and the Teleflex logo are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. © 2023 Teleflex Incorporated. All rights reserved.

"Rx only" is used within this labeling to communicate the following statement as presented in the FDA CFR : Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.

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