Inspected Dimensions:
Folded width: 5-1/2" (14 cm)
Folded length: 7-1/2" (19 cm)

ARROW®

Catheter Clamp and Fastener Product

Safety and Efficacy Considerations:

Do not use if package has been previously opened or damaged. **Warning: Prior to use read all package insert warnings, precautions, and instructions. Failure to do so may result in severe patient injury or death.**

Do not alter the catheter clamp or fastener during use or removal.

Procedure must be performed by trained personnel well versed in anatomical landmarks, safe technique, and potential complications.

Indications for Use:

The catheter clamp and fastener are indicated for use in securing catheter to patient or as a secondary suture site.

Contraindications:

None known.

Warnings and Precautions:*

- Warning: Sterile, Single use: Do not reuse, reprocess or resterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death.**
- Precaution: Catheter clamp and fastener must not be attached to catheter until spring-wire guide is removed.**

A Suggested Procedure: Use sterile technique.

- Follow manufacturer's instructions for placing catheter.
- After spring-wire guide has been removed and the necessary lines have been connected or locked, spread wings of rubber clamp and position on catheter as required to ensure proper tip location (refer to Fig. 1).


Fig. 1

- Snap rigid fastener onto catheter clamp (refer to Fig. 2).


Fig. 2

- Secure catheter to patient by suturing the catheter clamp and fastener together to the skin, using side wings to minimize the risk of catheter migration (refer to Fig. 3).


Fig. 3

- Continue procedure per manufacturer's instructions.

*If you have any questions, please contact Arrow International, Inc.

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* Om du har frågor kan du kontakta Arrow International, Inc.