## Peripheral Emergency Infusion Device - EID™

### Safety and Efficacy Considerations:

Do not use if package has been previously opened or damaged.

**Warning:** Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Do not alter EID™ or any kit/set component during insertion, use or removal.

Procedure must be performed by trained personnel well versed in anatomical landmarks, safe technique, and potential complications.

### Indications for Use:

The Arrow Peripheral Emergency Infusion Device (EID) permits access to the peripheral circulation for rapid volume infusion.

### Contraindications:

None known.

### Warnings and Precautions:\*

1. **Warning:** Sterile, Single use: Do not reuse, reprocess or resterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death.
2. **Warning:** Practitioners must be aware of complications associated with peripheral large bore catheters including vessel wall perforation, infiltration, catheter embolism, infection, thrombosis, inadvertent arterial puncture, nerve damage, hematoma, intravascular clotting, and hemorrhage.
3. **Warning:** Use appropriate size sheath for intended insertion site. Vessel rupture is possible if the diameter of the sheath is too large in comparison to the vessel diameter at the peripheral insertion site.
4. **Warning:** To lessen the risk of disconnects, only securely tightened Luer-Lock connections should be used with this device.
5. **Warning:** Do not apply excessive force in removing guide wire, dilator or sheath. If withdrawal cannot be easily accomplished, an x-ray should be obtained and further consultation requested.
6. **Warning:** Due to the risk of exposure to Human Immunodeficiency Virus (HIV) or other blood borne pathogens, health care workers should routinely use universal blood and body-fluid precautions in the care of all patients.
7. **Precaution:** Do not suture directly to the outside diameter of the sheath to minimize the risk of cutting or damaging the sheath or impeding sheath flow.
8. **Precaution:** Indwelling sheaths should be routinely inspected for desired flow rate, security of dressing. **Precaution:** Maintain the insertion site with regular meticulous redressing using aseptic technique.

### Rx only.

### A Suggested Procedure:

#### Use sterile technique.

1. Prep area of anticipated venipuncture.
2. Pre-assemble spring-wire guide tube to sheath/dilator assembly by placing male hub of tube into female hub of needle. Connection should be snug.
3. Use scalpel to nick skin at puncture site.
4. Do not remove clear guard from sheath/dilator assembly. Using guard tab, pull guard back toward hub to expose first 1/4 inch of sheath.
5. Puncture vessel using a continuous, controlled, slow, forward motion. **Precaution: Avoid transfixing both vessel walls. If both vessel walls are punctured, advancement of spring-wire guide could result in inadvertent subvessel placement.** Blood flashback in needle hub indicates positive placement in vein.  
The guard can be used to support the sheath during insertion. To remove guard continue to pull back. Guard will "peel off" sheath.  
For severely hypovolemic patients or patients in shock, it may be necessary to remove spring-wire guide assembly and attach a syringe to sheath assembly. Aspiration of syringe will hasten blood flashback. When placement of introducer in blood vessel is confirmed, syringe should be removed and spring-wire guide tube should be reconnected.
6. Stabilize position of sheath assembly and carefully advance spring-wire guide using actuating lever on spring-wire guide tube. Advance guide wire as far as required into blood vessel. **Precaution: If resistance is encountered while advancing spring-wire guide, do not force feed and do not retract spring-wire guide while in vessel. Withdraw entire unit and attempt new puncture.**
7. With wire in place rotate needle bevel 180° and advance assembly 1-2 mm further.
8. Maintain a firm grip on needle hub and advance sheath and dilator forward into vessel using slight twisting motion.
9. Hold sheath in place and remove dilator, needle, and spring-wire guide tube as a unit. Free blood flow indicates placement in vessel. **Warning: Do not leave tissue dilator in place as an indwelling catheter. Leaving tissue dilator in place puts patient at risk for possible vessel wall perforation. Warning: Do not reinsert needle into sheath/dilator assembly.**
10. Flush Luer-Lock infusion tubing and securely attach to sheath hub.
11. Secure sheath to patient. Suture wing is provided on hub. **Precaution: Do not suture directly to the outside diameter of the sheath to minimize the risk of cutting or damaging the sheath or impeding sheath flow.**
12. Cover puncture site with sterile dressing.
13. Record insertion procedure.

\*If you have any questions or would like additional reference information, please contact Arrow International, Inc.

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\*Om du har frågor eller önskar ytterligare referensinformation, kontakta Arrow International, Inc.