Peripheral Emergency Infusion Device - EID

For short term use (<30 days)

Rx only.

Indications for Use:

The Arrow Peripheral Emergency Infusion Device (EID) permits access to the peripheral circulation for rapid volume infusion.

Contraindications:

None known.

Clinical Benefits to be Expected:

The ability to access into the circulation and infuse large fluid volumes rapidly into a patient for treatment of shock or trauma, as examples.

⚠️ General Warnings and Precautions

Warnings:

1. Sterile, Single use: Do not reuse, reprocess or resterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death. Reprocessing of medical devices intended for single use only may result in degraded performance or a loss of functionality.
2. Read all package insert warnings, precautions and instructions prior to use. Failure to do so may result in severe patient injury or death.
3. Use appropriate size sheath for intended insertion site. Vessel rupture is possible if the diameter of the sheath is too large in comparison to the vessel diameter at the peripheral insertion site.
4. To lessen the risk of disconnects, only securely tightened Luer-Lock connections should be used with this device.
5. Do not apply excessive force in removing guidewire, dilator or sheath. If withdrawal cannot be easily accomplished, an x-ray should be obtained and further consultation requested.
6. Clinicians must be aware of complications/undesirable side-effects associated with peripheral large bore catheters including, but not limited to:

- vessel wall perforation
- infiltration
- catheter embolism
- infection
- thrombosis
- inadvertent arterial puncture
- nerve damage
- hematoma
- intravascular clotting
- hemorrhage

Precautions:

1. Do not alter the sheath or any other kit/set component during insertion, use, or removal.
2. Procedure must be performed by trained personnel well versed in anatomical landmarks, safe technique and potential complications.
3. Use standard precautions and follow institutional policies for all procedures including safe disposal of devices.
4. Do not suture directly to the outside diameter of the sheath to minimize the risk of cutting or damaging the sheath or impeding sheath flow.
5. Indwelling sheaths should be routinely inspected for desired flow rate, security of dressing,

6. Maintain the insertion site with regular meticulous redressing using aseptic technique.
7. Some disinfectants used at device insertion site contain solvents which can weaken the device material. Alcohol, acetone, and polyethylene glycol can weaken the structure of polyurethane materials. These agents may also weaken the adhesive bond between stabilization device and skin.
 - Do not use acetone on device surface.
 - Do not use alcohol to soak device surface or allow alcohol to dwell in a device lumen to restore patency or as an infection prevention measure.
 - Do not use polyethylene glycol containing ointments at insertion site.
 - Take care when infusing drugs with a high concentration of alcohol.
 - Allow insertion site to dry completely prior to applying dressing.

A Suggested Procedure: Use Sterile Technique.

⚠️ Precaution: Place patient in slight Trendelenburg position as tolerated to reduce the risk of air embolism. If femoral approach is used, place patient in supine position.

1. Prep area of anticipated venipuncture.
2. Pre-assemble guidewire tube to sheath/dilator assembly by placing male hub of tube into female hub of needle. Connection should be snug.
3. Use scalpel to nick skin at puncture site.
4. Do not remove clear guard from sheath/dilator assembly. Using guard tab, pull guard back toward hub to expose first 1/4 inch of sheath.
5. Puncture vessel using a continuous, controlled, slow, forward motion.

⚠️ Precaution: Avoid transfixing both vessel walls. If both vessel walls are punctured, advancement of guidewire could result in inadvertent subvessel placement.

Blood flashback in needle hub indicates positive placement in vein.

The guard can be used to support the sheath during insertion. To remove guard continue to pull back. Guard will "peel off" sheath.

For severely hypovolemic patients or patients in shock, it may be necessary to remove guidewire assembly and attach a syringe to sheath assembly. Aspiration of syringe will hasten blood flashback. When placement of introducer in blood vessel is confirmed, syringe should be removed and guidewire tube should be reconnected.

6. Stabilize position of sheath assembly and carefully advance guidewire using actuating lever on guidewire tube. Advance guidewire as far as required into blood vessel.

⚠️ Precaution: If resistance is encountered while advancing guidewire, do not force feed and do not retract guidewire while in vessel. Withdraw entire unit and attempt new puncture.

7. With wire in place rotate needle bevel 180° and advance assembly 1-2 mm further.
8. Maintain a firm grip on needle hub and advance sheath and dilator forward into vessel using slight twisting motion.
9. Hold sheath in place and remove dilator, needle, and guidewire tube as a unit. Free blood flow indicates placement in vessel.

⚠️ Warning: Do not leave tissue dilator in place as an indwelling catheter. Leaving tissue dilator in place puts patient at risk for possible vessel wall perforation.

⚠️ Warning: Do not reinsert needle into sheath/dilator assembly.

10. Flush Luer-Lock infusion tubing and securely attach to sheath hub.
 11. Secure sheath to patient. Suture wing is provided on hub.
- ⚠️ Precaution: Do not suture directly to the outside diameter of the sheath to minimize the risk of cutting or damaging the sheath or impeding sheath flow.**
12. Cover puncture site with sterile dressing.
 13. Record insertion procedure.

For reference literature concerning patient assessment, clinician education, insertion technique, and potential complications associated with this procedure, consult standard textbooks, medical literature, and Arrow International LLC website: www.teleflex.com

A pdf copy of this IFU is located at www.teleflex.com/IFU

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority. The contacts of national competent authorities (Vigilance Contact Points) and further information can be found on the following European Commission website: https://ec.europa.eu/growth/sectors/medical-devices/contacts_en

 Symbol Glossary: Symbols are in compliance with ISO 15223-1. Some symbols may not apply to this product. Refer to product labeling for symbols that apply specifically to this product.								
								
Caution	Medical device	Consult instructions for use	Do not reuse	Do not resterilize	Sterilized by ethylene oxide	Single sterile barrier system	Keep away from sunlight	Keep dry
								
Do not use if package is damaged	Not made with natural rubber latex	Catalogue number	Lot number	Use by	Manufacturer	Date of manufacture	Importer	

Arrow, the Arrow logo, Teleflex and Teleflex logo are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. © 2021 Teleflex Incorporated. All rights reserved.

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