# Rapid Infusion Catheter (RIC®) Exchange Product

For short term use (<30 days)

Rx only

## Indications for Use:

The Arrow Rapid Infusion Catheter (RIC®) exchange product permits exchange of a 20 Ga. (or larger) peripheral intravenous catheter to a large diameter vascular access device for rapid volume infusion.

## Contraindications:

None known.

## Clinical Benefits to be Expected:

The ability to provide rapid access to the vascular system in the trauma or emergency situation.

## ⚠️ General Warnings and Precautions

### Warnings:

1. Sterile, Single use: Do not reuse, reprocess or resterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death. Reprocessing of medical devices intended for single use only may result in degraded performance or a loss of functionality.
2. Read all package insert warnings, precautions and instructions prior to use. Failure to do so may result in severe patient injury or death.
3. Use appropriate size device for intended exchange site. Vessel rupture is possible if the diameter of the device is too large in comparison to the vessel diameter at the peripheral exchange site.
4. To lessen the risk of disconnects, only securely tightened Luer-Lock connections should be used with this device.
5. Do not apply excessive force in removing guidewire, dilator or vascular access device. If withdrawal cannot be easily accomplished, an x-ray should be obtained and further consultation requested.
6. Clinicians must be aware of complications/undesirable side effects associated with peripheral large bore catheters including, but not limited to:
  - vessel wall perforation
  - nerve damage
  - infiltration
  - hematoma
  - catheter embolism
  - intravascular clotting
  - infection
  - hemorrhage
  - thrombosis
  - inadvertent arterial puncture

### Precautions:

1. Do not alter any kit/set component during insertion, use, or removal.
2. Procedure must be performed by trained personnel well versed in anatomical landmarks, safe technique and potential complications.
3. Use standard precautions and follow institutional policies for all procedures including safe disposal of devices.
4. Do not suture directly to the outside diameter of the vascular access device to minimize the risk of cutting or damaging the device or impeding flow.

5. Vascular access devices should be routinely inspected for desired flow rate, security of dressing, correct position and for proper Luer-Lock connection.
6. Some disinfectants used at device insertion site contain solvents which can weaken the device material. Alcohol, acetone, and polyethylene glycol can weaken the structure of polyurethane materials. These agents may also weaken the adhesive bond between stabilization device and skin.
  - Do not use acetone on device surface.
  - Do not use alcohol to soak device surface or allow alcohol to dwell in a device lumen to restore patency or as an infection prevention measure.
  - Do not use polyethylene glycol containing ointments at insertion site.
  - Take care when infusing drugs with a high concentration of alcohol.
  - Allow insertion site to dry completely prior to applying dressing.
7. Maintain the insertion site with regular meticulous redressing using aseptic technique.

**Kits/Sets may not contain all accessory components detailed in these instructions for use. Become familiar with instructions for individual component(s) before beginning the procedure.**

## A Suggested Procedure: Use Sterile Technique.

1. Position patient as appropriate for insertion site.
  2. Prepare clean skin with an appropriate antiseptic agent.
  3. Drape puncture site.
  4. Administer local anesthetic per institutional policies and procedures.
  5. Disconnect tubing from catheter and insert guidewire the required distance into catheter.
- ⚠️ **Precaution: Maintain firm grip on guidewire at all times.**
- ⚠️ **Precaution: If resistance is encountered while advancing guidewire, hold catheter in place and carefully withdraw guidewire. Attempt to reinsert. If resistance is encountered again, abort procedure.**
6. With guidewire in place within the vein, remove indwelling peripheral catheter over guidewire.
  7. Enlarge cutaneous puncture site with cutting edge of scalpel positioned away from guidewire.
- ⚠️ **Precaution: Do not cut guidewire.**
8. Remove plastic guard from vascular access device/dilator assembly.
  9. Thread tapered tip of dilator over guidewire. Grasping near skin, advance dilator and vascular access device into vessel with a slight twisting motion.
  10. Advance device over dilator into vessel, again grasping near skin and using slight twisting motion.
  11. Holding device in place, remove guidewire and dilator as a unit.
- ⚠️ **Precaution: Do not leave vessel dilator in place as an indwelling catheter to minimize the risk of possible vessel wall perforation.**
12. Attach desired Luer-Lock administration tubing to device hub. For best results large diameter fluid administration tubing should be used.

**⚠️ Precaution: To minimize the risk of disconnects use only Luer-Lock connecting tubing.**

13. Secure vascular access device to patient using suture wing and/or anchor with a purse string suture around the device suture ring.

**⚠️ Precaution: Do not suture directly to the outside diameter of the vascular access device to minimize the risk of cutting or damaging the device or impeding flow.**

14. Dress puncture site per hospital protocol.

15. Record the exchange procedure on the patient's chart.

For reference literature concerning patient assessment, clinician education, insertion technique, and potential complications associated with this procedure, consult standard textbooks, medical literature, and Arrow International LLC website: [www.teleflex.com](http://www.teleflex.com)

A pdf copy of this IFU is located at [www.teleflex.com/IFU](http://www.teleflex.com/IFU)

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority. The contacts of national competent authorities (Vigilance Contact Points) and further information can be found on the following European Commission website: [https://ec.europa.eu/growth/sectors/medical-devices/contacts\\_en](https://ec.europa.eu/growth/sectors/medical-devices/contacts_en)

**en**

**Symbol Glossary: Symbols are in compliance with ISO 15223-1.**  
Some symbols may not apply to this product. Refer to product labeling for symbols that apply specifically to this product.

								
Caution	Medical device	Consult instructions for use	Do not reuse	Do not resterilize	Sterilized by ethylene oxide	Single sterile barrier system	Keep away from sunlight	Keep dry
								
Do not use if package is damaged	Not made with natural rubber latex	Catalogue number	Lot number	Use by	Manufacturer	Date of manufacture		

*Arrow, the Arrow logo, RIC, Teleflex and Teleflex logo are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. © 2024 Teleflex Incorporated. All rights reserved.*

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