# Introducer Needle Products/ Catheter Over Introducer Needle Products

Rx only

## Intended Purpose:

These devices are intended to provide access to the vascular system for the introduction of a guidewire.

## Indications for Use:

The Introducer Needle allows access to the vascular system for the introduction of a guidewire to facilitate catheter placement.

The Introducer Catheter over Needle is intended to puncture the skin to gain access to the vascular system and, after removal of the introducer needle, allows passage of a spring wire guide into the vessel.

## Patient Target Group:

Intended to be used in patients with anatomy suitable for use with the device.

## Contraindications:

None known.

## Clinical Benefits to be Expected:

Introducer Needle: The ability to provide access to a vessel with a single puncture site for introduction of a guidewire and catheter.

Catheter Over Needle: The ability to puncture the skin to gain access to the vascular system and, after removal of the introducer needle, to allow passage of a spring wire guide into the vessel.

## Performance Characteristics:

Introducer Needle: Access to a vessel with a single puncture site for introduction of a guidewire.

Catheter over Introducer Needle: Allows passage of a guidewire into the vessel for catheter introduction.

## Contains Hazardous Substance:

Components manufactured using Stainless Steel can contain > 0.1% weight by weight of Cobalt (CAS # 7440-48-4) which is considered a category 1B CMR (Carcinogenic, mutagenic or toxic to reproduction) substance. The amount of Cobalt in the Stainless Steel components has been evaluated and considering the intended purpose and toxicological profile of the devices there is no biological safety risk to patients when using the devices as instructed within this IFU.

## General Warnings and Precautions

### Warnings:

1. Sterile, Single use: Do not reuse, reprocess or resterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death. Reprocessing of medical devices intended for single use may result in degraded performance or a loss of functionality.
2. Read all package insert warnings, precautions and instructions prior to use. Failure to do so may result in severe patient injury or death.

3. Air embolism can occur if air is allowed to enter a central venous access device or vein. Do not leave open needles in central venous puncture site. Use only securely tightened Luer-Lock connections with any central venous access device to guard against inadvertent disconnection.

4. Clinicians must be aware of complications/undesirable side-effects associated with introducer needles including, but not limited to:

- Particle/Foreign Embolism
- Local Tissue Inflammation
- Bleeding/Hematoma/ Hemorrhage
- Vessel perforation
- Vessel trauma
- Bloodstream Infection
- External Bleeding
- Local Infection
- Arterial puncture
- Phlebitis
- Pneumothorax
- Hemothorax
- Thrombosis
- Accidental needlesticks
- Exposure to contaminated blood

## Precautions:

1. Procedure must be performed by trained personnel well versed in anatomical landmarks, safe technique and potential complications.
2. Use standard precautions and follow institutional policies for all procedures including safe disposal of devices.
3. Some disinfectants used at catheter insertion site contain solvents which can weaken the catheter material. Alcohol, acetone, and polyethylene glycol can weaken the structure of polyurethane materials. These agents may also weaken the adhesive bond between catheter stabilization device and skin.
  - Do not use acetone on catheter surface.
  - Do not use alcohol to soak catheter surface or allow alcohol to dwell in a catheter lumen to restore catheter patency or as an infection prevention measure.
  - Do not use polyethylene glycol containing ointments at insertion site.
  - Take care when infusing drugs with a high concentration of alcohol.
  - Allow insertion site to dry completely prior to procedure and before applying dressing.
  - Do not allow kit components to come into contact with alcohol.
4. If the package is damaged or unintentionally opened before use do not use the device. Dispose of the device.
5. Storage conditions for these devices require that they are kept dry and out of direct sunlight.

## A Suggested Procedure: Use sterile technique.

### Prep Puncture Site:

1. Prep the puncture site per institutional policies and procedures.

- Prepare and clean skin with an appropriate antiseptic agent and allow to dry.

## Gain Initial Venous Access:

### Echogenic Needle (where provided):

An echogenic needle is used to allow access to the vascular system for the introduction of a guidewire to facilitate catheter placement. The needle tip is enhanced for approximately 1 cm for clinician to identify exact needle tip location when puncturing the vessel under ultrasound.

### Protected Needle/Safety Needle (where provided):

A protected needle/safety needle should be used in accordance with manufacturer's instructions for use.

**⚠ Warning:** Do not leave open needles or catheters in central venous puncture site. Air embolism can occur if air is allowed to enter a central venous access device or vein.

**⚠ Precaution:** Do not reinsert needle into introducer catheter (where provided) to reduce risk of catheter embolus.

- Locate the vein as specified in catheter associated instructions for use and aspirate.

## Verify Venous Access:

- Verify correct venous access as specified in catheter associated instructions for use. For catheter over needle use, remove the needle leaving the catheter in place.

**⚠ Warning:** Pulsatile flow is usually an indicator of inadvertent arterial puncture.

**⚠ Precaution:** Do not rely on blood aspirate color to indicate venous access.

## Insert Guidewire:

- Continue with guidewire introduction and catheter placement procedure as specified in catheter associated instructions for use.

**⚠ Warning:** Do not withdraw guidewire against needle bevel to reduce risk of possible severing or damaging of guidewire.

For reference literature concerning patient assessment, clinician education, insertion technique, and potential complications associated with this procedure, consult standard textbooks, medical literature, and Arrow International LLC website: [www.teleflex.com](http://www.teleflex.com)

A pdf copy of this IFU is located at [www.teleflex.com/IFU](http://www.teleflex.com/IFU)

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority. The contacts of national competent authorities (Vigilance Contact Points) and further information can be found on the following European Commission website: [https://ec.europa.eu/growth/sectors/medical-devices/contacts\\_en](https://ec.europa.eu/growth/sectors/medical-devices/contacts_en)

en

Symbol Glossary: Symbols are in compliance with ISO 15223-1.

Some symbols may not apply to this product. Refer to product labeling for symbols that apply specifically to this product.

Caution	Medical device	Consult instructions for use	Contains a medicinal substance	Contains hazardous substances	Do not reuse	Do not resterilize	Sterilized by ethylene oxide
Single sterile barrier system with protective packaging inside	Single sterile barrier system	Keep away from sunlight	Keep dry	Do not use if package is damaged	Not made with natural rubber latex	Catalogue number	Lot number
				<i>Arrow, the Arrow logo, Teleflex and the Teleflex logo are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. © 2023 Teleflex Incorporated. All rights reserved.</i>			
Use by	Manufacturer	Date of manufacture	Importer	<i>"Rx only" is used within this labeling to communicate the following statement as presented in the FDA CFR: Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.</i>			

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