

CG⁺™ ARROW® PICC Pressure Injectable Information

MAXIMUM injector pressure setting: 300 psi

Cross Section	Catheter Size	Lumen	Priming Volume (mL) ¹	Gravity Flow Rate (mL/hr) ²	Pump Flow Rate (mL/hr) ³	MAX Indicated Pressure Injection Flow Rate (11.8 cP) ^{4,5}	MAX Catheter Pressure During MAX Flow Conditions ^{5,6}	Minimum Static Burst Pressure ^{5,7}
	4.5 Fr. 40 cm 1-L	Distal (D)	0.49	1510	6350	5 mL/sec	210 psi	411 psi
	4.5 Fr. 45 cm 1-L	Distal (D)	0.50	1360	6290	5 mL/sec	223 psi	460 psi
	4.5 Fr. 50 cm 1-L	Distal (D)	0.57	1210	5870	5 mL/sec	248 psi	382 psi
	4.5 Fr. 55 cm 1-L	Distal (D)	0.59	1100	5130	5 mL/sec	257 psi	423 psi
	5.5 Fr. 40 cm 2-L	Distal (D)	0.41	500	3800	5 mL/sec	206 psi	352 psi
		Proximal (P)	0.43	520	3770	5 mL/sec	207 psi	352 psi
	5.5 Fr. 45 cm 2-L	Distal (D)	0.43	520	3230	5 mL/sec	218 psi	402 psi
		Proximal (P)	0.45	560	3410	5 mL/sec	215 psi	414 psi
	5.5 Fr. 50 cm 2-L	Distal (D)	0.45	400	3180	5 mL/sec	220 psi	362 psi
		Proximal (P)	0.47	410	3150	5 mL/sec	220 psi	357 psi
	5.5 Fr. 55 cm 2-L	Distal (D)	0.49	360	3020	5 mL/sec	230 psi	346 psi
		Proximal (P)	0.50	370	2970	5 mL/sec	227 psi	323 psi
	6 Fr. 40 cm 3-L	Distal (D)	0.50	1370	6970	6 mL/sec	187 psi	357 psi
		Medial (M)	0.39	490	2810	not rated	not rated	not rated
		Proximal (P)	0.41	510	2880	not rated	not rated	not rated
	6 Fr. 45 cm 3-L	Distal (D)	0.54	1190	6410	6 mL/sec	187 psi	353 psi
		Medial (M)	0.41	460	2610	not rated	not rated	not rated
		Proximal (P)	0.42	450	2620	not rated	not rated	not rated
	6 Fr. 50 cm 3-L	Distal (D)	0.57	1120	5920	6 mL/sec	196 psi	351 psi
		Medial (M)	0.43	410	2370	not rated	not rated	not rated
		Proximal (P)	0.44	420	2390	not rated	not rated	not rated
	6 Fr. 55 cm 3-L	Distal (D)	0.61	1060	5360	6 mL/sec	207 psi	349 psi
		Medial (M)	0.45	370	2090	not rated	not rated	not rated
		Proximal (P)	0.47	390	2160	not rated	not rated	not rated

The approximate flow capability for the 4.5 Fr. 1-L and 5.5 Fr. 2-L pressure injectable lumens using an injectate of 4.7 centipoise (cP) is 8 mL/sec. The approximate flow capability for the 6 Fr. 3-L pressure injectable lumen using an injectate of 4.7 cP is 9 mL/sec. ¹Priming volumes are approximate and done without injection caps. ²Flow rates were determined using room temperature water, 100 cm head height and represent approximate flow capabilities. ³Pump flow rates are determined at maximum pump pressure of 10 psig and represent approximate flow capabilities. ⁴Pressure injection flow rates are determined at the injector pressure setting of 300 psi maximum using media of 11.8 centipoise viscosity, with 152 cm pressure tubing. ⁵Based on testing conducted using media at a viscosity of 11.8 cP measured at body temperature (37 °C). ⁶Internal catheter pressure during pressure injection at maximum flow conditions with injector safety cut-off at 300 psi. ⁷Minimum static burst pressure is the failure point of catheter when occluded.

⚠ Pressure Injection General Warnings and Cautions:

Warnings:

1. Assess each patient for appropriateness of a pressure injection procedure.
2. Pressure injection procedures must be performed by trained personnel well versed in safe technique and potential complications.
3. Discontinue pressure injections at first sign of infiltration/extravasation. Follow hospital/institutional protocol for appropriate medical intervention.

Cautions:

1. Do not exceed ten (10) injections or the maximum pressure of 300 psi on power injector equipment to reduce the risk of catheter failure and/or tip displacement.
2. Pressure limit settings on power injector equipment may not prevent over pressurizing an occluded or partially occluded catheter.
3. Use appropriate administration set tubing between catheter and power injector equipment to reduce the risk of catheter failure.

Pressure Injection Instructions:

Use sterile technique.




⚠ Warning: Use an appropriate method to confirm catheter tip position prior to each pressure injection per hospital/institutional policy.

1. Remove injection cap from the extension line of catheter lumen to be injected.
2. Check for patency through intended pressure injectable lumen:
 - Attach 10 mL syringe, or larger, filled with sterile normal saline.
 - Flush catheter.
 - Aspirate catheter until approximately 3 mL of blood enters syringe freely.
 - Vigorously flush catheter.

⚠ Warning: Ensure patency of intended pressure injectable lumen of catheter prior to pressure injection to reduce the risk of catheter failure and/or patient complications.

3. Detach syringe.
4. Attach pressure injection administration set tubing to appropriate extension line of catheter according to manufacturer's recommendations.
 - ⚠ Warning: Use only lumen(s) labeled "Pressure Injectable" for pressure injection to reduce the risk of catheter failure and/or patient complications.
 - ⚠ Caution: Do not exceed ten (10) injections or catheter's maximum recommended flow rate located on product labeling to reduce the risk of catheter failure and/or tip displacement.
5. Inject contrast media in accordance with hospital/institutional protocol.
 - ⚠ Caution: Warm contrast media to body temperature prior to pressure injection to reduce the risk of catheter failure.
 - ⚠ Caution: Follow the contrast media manufacturer's specified instructions for use, contraindications, warnings and precautions.
6. Disconnect catheter from power injector equipment.
7. Flush catheter using a 10 mL syringe, or larger, filled with sterile normal saline.
8. Disconnect syringe and replace with sterile injection cap on catheter extension line.

 See Instructions for Use for additional information.

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