VPS Rhythm™ Device ECG Accessories

Only for use with the VPS Rhythm™ Device

Rx only

Intended Use:

The VPS Rhythm™ Device is intended to provide real time tip location information of a central venous catheter by utilization of ECG to observe P-wave changes as the tip approaches the right atrium of the heart via the Superior Vena Cava (SVC).

The VPS Rhythm Device is appropriate for final catheter tip placement in any catheter that is intended to be placed in the lower 1/3 of the SVC near the cavo-atrial junction (CAJ). Catheters that are traditionally considered CVCs that are intended to be placed in the right atrium, or any location other than the lower 1/3 SVC (e.g. chronic hemodialysis catheters placed in the right atrium), must be verified using another method.

Indications for Use for the VPS Rhythm Device:

The VPS Rhythm Device is indicated for the positioning of central venous catheters including PICCs. It provides catheter tip location information by using the patient's cardiac electrical activity. The VPS Rhythm Device is indicated for use as an alternative method to chest X-ray or fluoroscopy for confirmation of central venous catheter tip placement in adult patients. The TipTracker Technology is an optional accessory for use with the VPS Rhythm Device, indicated for visual navigation of a peripherally-inserted central catheter (PICC) as it is inserted through the vasculature. The TipTracker technology is used for catheter tip navigation purposes only; it is not used to determine final catheter tip placement.

NOTE: In general, devices that utilize ECG technique to observe P-wave are limited, but not contraindicated, for patients where cardiac rhythms may change presentation of the P-wave; including:

- Atrial fibrillation
- Pacemaker-driven rhythm
- Atrial flutter
- Chronic obstructive pulmonary disease (COPD)
- Severe tachycardia

Such patients are easily identified prior to central catheter insertion. In these specific cases, use of an additional confirmation method is necessary to confirm catheter tip location.

Device Description:

The VPS Rhythm Device is designed to guide the clinician in central venous catheter tip placement to an area within the lower 1/3 of the superior vena cava (SVC), near the cavo-atrial junction (CAJ). The VPS Rhythm Device can be used with various types and brands of central venous catheters, including PICCs, CVCs, Ports and Acute Hemodialysis catheters.

The system uses two ECG leads placed on the patient's chest and utilizes a third (intravascular) ECG lead from a wire inside the catheter or from a column of normal saline through the catheter lumen. As the catheter approaches the heart, electrocardiograph (ECG) waveforms are used to determine final catheter tip location.

For the use of the VPS Rhythm Device, two sterile, single use components are used: a sterile sleeve for the Remote Control and cable and an ECG Clip Cable. The sterile sleeve provides a means for using the VPS Rhythm Device Remote Control to operate the VPS Rhythm monitor within the sterile field. The ECG Clip Cable is used to connect the VPS Rhythm Device to a wire inside the catheter lumen. The ECG Clip Cable may optionally be connected to an Arrow-Johans™ adapter (where provided) to connect the VPS Rhythm Device to a column of normal saline through the catheter lumen.

Contraindications:

There are no contraindications associated with the use of the ECG accessories. Consult the catheter manufacturer's Instructions for Use for possible catheter contraindications.

Clinical Benefits to be Expected:

The VPS Rhythm Device's clinical benefits are the following:

Makes central venous catheter insertion safer by reducing procedure-related complications, Improves efficiency and speed when inserting central venous catheters, increasing first-attempt success rate compared to "blind" insertions without navigation and/or tip confirmation, provides safe and reliable real-time guidance during insertion of PICCs (when TipTracker is being used) to reduce catheter misplacement, provides confirmation of final catheter tip location to reduce catheter malposition, reducing treatment delays due

to reliance on X-rays confirmation of catheter tip location, and reduces radiation exposure due to use of X-rays to confirm catheter tip location.

⚠ General Warnings and Precautions

Warnings:

1. **Sterile, Single use: Do not reuse, reprocess or resterilize.**
Reuse of device creates a potential risk of serious injury and/or infection which may lead to death. Reprocessing of medical devices intended for single use only may result in degraded performance or a loss of functionality.
2. Read all package insert warnings, precautions and instructions prior to use. Failure to do so may result in severe patient injury or death.
3. Refer to the VPS Rhythm Device Operator's Manual and this IFU for full instructions, indications, contraindications, warnings and precautions. Failure to do so may result in severe patient injury or death.

Precautions:

1. Do not alter any kit/set component during insertion, use or removal (except as instructed).
2. Procedure must be performed by trained personnel well versed in anatomical landmarks, safe technique and potential complications.
3. Use standard precautions and follow institutional policies for all procedures including safe disposal of devices.

Preparation and Use of the ECG Accessories with the VPS Rhythm Device:

Follow institutional policy and procedure and manufacturer's guidelines for catheter placement. The catheter should be inserted, manipulated, and removed by a qualified, licensed physician or other qualified health care professional, familiar with the use of the VPS Rhythm Device.

Non-Sterile Setup:

1. Turn on the VPS Rhythm Device and follow full instructions in the VPS Rhythm Device Operator's Manual for entering the patient ID and entering optional notes.
2. Connect black ECG Cable connector to the black socket and the blue Remote Control connector to the blue socket on the back of the monitor.
3. Place ECG electrodes.
 - a. Attach ECG snap leads to corresponding color of the ECG Cable hub and attach the electrodes to the snap leads.
 - b. Prepare the skin according to institutional policy and procedure to ensure good adherence and electrical contact. Ensure the electrodes are not over any bony prominence.

AHA ECG Cable Connections and Placement

- White to right arm (RA)
- Black to left arm (LA)
- Red to lower left chest or upper left leg (LL)

IEC ECG Cable Connections and Placement

- Red to right arm (R)
- Yellow to left arm (L)
- Green to lower left chest or upper left leg (F)

4. Follow VPS Rhythm Operator's Manual to obtain an External ECG baseline. Confirm patient is in sinus rhythm.
5. If using PICC mode, once baseline external ECG is saved, a window will appear to allow user to enter external measurement and trimmed length. Using landmark technique, obtain external measurement for PICC length; enter into external measurement window. Use Remote Control or touch screen to enter external measurement and trimmed length.

Add 2 cm to external measured length and enter as trimmed length. Use this length for trimming, if trimming is required.

Sterile Procedure:

6. Prepare sterile field per institutional policy and procedure and manufacturer's guidelines.
7. Connect ECG Clip Cable to Remote Control and cable.
 - a. The sterile sleeve is similar to an ultrasound probe sleeve and is folded showing an arrow at one end and instructional images at the other. Starting at the instructional end, locate the hand illustration inside the center fold of the sleeve and carefully place your hand all the way inside.
 - b. Maintaining sterile technique, firmly grasp the tip of the Remote cable plug through the end of the sterile sleeve where the arrow is located and unroll the sterile sleeve to cover the entire length of the Remote Control and cable.
 - c. Pull the plastic sleeve tightly over the tip of Remote Control cable plug creating tension against tip of plug.
 - d. Connect the Remote Control plug to the ECG Clip Cable jack, carefully piercing through the sterile sleeve. **Ensure sterile field is maintained.**
8. Prepare catheter per institutional policy and procedure and catheter manufacturer's Instructions for Use.
 - a. Ensure monitor is in Intravascular mode.

Using Valved Catheters

If using a valved catheter special considerations must be made for correct use of the VPS Rhythm Device for catheter tip placement. The ECG signal used for final catheter tip placement is impeded if the catheter's integral valve is closed.

- *Catheters with a valve at the distal end: Use the saline column technique with Johans adapter to insert these catheters. While advancing catheter towards the SVC, periodically flush the catheter through the distal end of the Johans adapter to open the PICC's integral valve.*
- *Catheters with a valve integrated into the catheter's luer hub:

 - *The saline column technique with Johans adapter can be used: While advancing the catheter towards the SVC, periodically flush the catheter through the distal end of the Johans Adapter to open the catheter's integral valve.*
 - *If the catheter has a metal stylet/wire where the position of the wire can be secured relative to the catheter hub, the technique below may be used. When using this technique the metal stylet/wire passes through the valve, keeping it open throughout the procedure.**

NOTE: Always maintain a complete saline column inside the catheter.

NOTE: To maintain a clear waveform, it is recommended to keep the patient as still as possible.

9. Insert catheter.

For catheters with a metal stylet/wire, when the position of the stylet/wire can be secured relative to the catheter hub:

- a. Flush all catheter lumens with sterile saline and insert catheter wire into lumen so that the tip of the wire is within the catheter.

- b. Connect ECG Clip Cable to catheter wire near the catheter's luer hub by squeezing the alligator clip and placing the clip over the metal portion of the catheter wire/stylet.
- c. Use Modified Seldinger technique, inserting the catheter and wire into the vessel as a unit. Place the catheter per institutional policy and procedure. Thread the catheter slowly.
- d. As the catheter approaches the SVC, observe changes in the ECG waveform to place the line at the patient's maximum P-wave, the proper tip location. Refer to the VPS Rhythm Device Operator's Manual for instruction on catheter tip placement and tip confirmation, including using Intravascular ECG.

For catheters with a non-metal stylet: Saline column technique using Johans adapter

- a. While maintaining sterile technique, connect the Johans Adapter to the catheter via the luer lock.
- b. Remove the rubber protective cap from the Johans Adapter.
- c. Connect the ECG Clip Cable to the catheter by squeezing the ends of the alligator clip and connecting it to the exposed metal on the Johans Adapter.
- d. Flush the Johans Adapter and all catheter lumens with sterile saline. Replace the white cap on the distal Johans Adapter luer connector.
- e. Place the central line using the Seldinger technique per institutional policy and procedure and manufacturer's guidelines.
- f. As the catheter approaches the SVC, observe changes in the ECG waveform to place the line at the patient's maximum P-wave, the proper tip location. Refer to the VPS Rhythm Device Operator's Manual for instructions on catheter tip placement and tip confirmation, including using Intravascular ECG.

Complete Procedure:

10. Remove ECG Clip Cable (or Johans adapter) from catheter. Complete procedure per institutional policy and procedure and catheter manufacturer's Instructions for Use.

For reference literature concerning patient assessment, clinician education, insertion technique, and potential complications associated with this procedure, consult standard textbooks, medical literature, and Arrow International LLC website: www.teleflex.com

A pdf copy of this IFU is located at www.teleflex.com/IFU

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority. The contacts of national competent authorities (Vigilance Contact Points) and further information can be found on the following European Commission website: https://ec.europa.eu/growth/sectors/medical-devices/contacts_en

en

Symbol Glossary: Symbols are in compliance with ISO 15223-1.

Some symbols may not apply to this product. Refer to product labeling for symbols that apply specifically to this product.

								
Caution	Medical device	Consult instructions for use	Do not reuse	Do not resterilize	Sterilized by ethylene oxide	Single sterile barrier system	Keep away from sunlight	Keep dry
								
Do not use if package is damaged	Not made with natural rubber latex	Catalogue number	Lot number	Use by	Manufacturer	Date of manufacture	Importer	

Arrow, the Arrow Logo, Arrow-Johans, Teleflex, the Teleflex Logo, TipTracker and VPS Rhythm are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. © 2024 Teleflex Incorporated. All rights reserved.

"Rx only" is used within this labeling to communicate the following statement as presented in the FDA CFR : Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.

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