# ECG (Electrocardiogram) Extension Cable

## Rx only

### Intended Purpose/Indications for Use:

The ECG Extension Cable permits verification of central venous catheter placement by right atrial electrophysiology.

### Patient Target Group:

Intended for use when an ECG adapter is used for central venous catheter placement.

### Contraindications:

None known.

### Clinical Benefits to be Expected:

The ECG Extension Cable provides confirmation of central venous catheter placement by right atrial electrophysiology.

### Performance Characteristics:

This device provides a method to confirm placement of central venous catheter by right atrial electrophysiology.

### ⚠️ General Warnings and Precautions

#### Warnings:

1. Sterile, Single use: Do not reuse, reprocess or resterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death. Reprocessing of medical devices intended for single use only may result in degraded performance or a loss of functionality.
2. Read all package insert warnings, precautions and instructions prior to use. Failure to do so may result in severe patient injury or death.
3. Strict electrical isolation and grounding techniques must be used for monitoring equipment.
4. Do not leave cable attached to guidewire when cable not attached to monitor. An electrical shock hazard is created due to a direct electrical path to the patient. Always detach from guidewire first when disconnecting ECG extension cable.

#### Precautions:

1. Do not alter the cable during use.
2. Procedure must be performed by trained personnel well versed in anatomical landmarks, safe technique and potential complications.
3. Use standard precautions and follow institutional policies for all procedures including safe disposal of devices.
4. If the package is damaged or unintentionally opened before use do not use the device. Dispose of the device.
5. Storage conditions for these devices require that they are kept dry and out of direct sunlight.

### A Suggested Procedure: Use sterile technique.

#### 1. a) for Arrow catheters 20 and 30 cm:

Place the catheter into vessel using over-the-wire technique per manufacturer's instructions. When using a marked guidewire, refer to Table 1 to determine when the wire is just outside the catheter tip based on Arrow® CVC catheter lengths.

20 cm catheter	approximately 35 cm wide mark	information related to 60 cm ECG guidewire
30 cm catheter	approximately 45 cm wide mark	information related to 68 cm ECG guidewire

Table 1


Figure 1

#### 1. b) for Arrow catheters 16 cm:

Place the catheter into vessel using over-the-wire technique per manufacturer's instructions. When 32 cm mark on guidewire exits 2 cm from catheter hub, the wire "J" is just outside the catheter tip.


Figure 2

#### 1. c) for non-Arrow catheters:

1. Follow instructions for use of the manufacturer of the catheter and guidewire.
2. Connect cable to monitor.
3. Connect alligator clamp of ECG extension cable to guidewire at referenced mark as described above.

#### ⚠️ Precaution: Do not allow guidewire to move within the catheter.

4. Advance catheter and guidewire as a unit until P wave of ECG trace begins to progressively get larger. At this point, the catheter tip should be at the cavo-atrial junction. If the tip is advanced further into the right atrium, the P wave will either be as large as the R wave or become biphasic.  
For catheters 16 cm or shorter the catheter may not reach the right atrium, and this technique can only be used to determine if the catheter is facing the proper direction. Attach extension cable as described above and advance catheter and guidewire as a unit. When the P wave becomes higher than normal the catheter is facing the proper direction.
5. Withdraw catheter and guidewire as a unit 2-3 cm from cavo-atrial junction. At this position the tip should be within the superior vena cava.
6. Holding catheter in place remove guidewire.
7. Continue procedure per manufacturer's instructions.

For reference literature concerning patient assessment, clinician education, insertion technique, and potential complications associated with this procedure, consult standard textbooks, medical literature, and Arrow International LLC website: [www.tелефex.com](http://www.tелефex.com)

A pdf copy of this IFU is located at [www.teleflex.com/IFU](http://www.teleflex.com/IFU)

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority. The contacts of national competent authorities (Vigilance Contact Points) and further information can be found on the following European Commission website: [https://ec.europa.eu/growth/sectors/medical-devices/contacts\\_en](https://ec.europa.eu/growth/sectors/medical-devices/contacts_en)

en Symbol Glossary: Symbols are in compliance with ISO 15223-1								
								
Caution	Medical device	Consult instructions for use	Do not reuse	Do not resterilize	Sterilized by ethylene oxide	Single sterile barrier system	Keep away from sunlight	Keep dry
								
Do not use if package is damaged	Not made with natural rubber latex	Catalogue number	Lot number	Use by	Manufacturer	Date of manufacture	Importer	

*Arrow, the Arrow logo, Teleflex and the Teleflex logo are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. © 2024 Teleflex Incorporated. All rights reserved.*

*"Rx only" is used within this labeling to communicate the following statement as presented in the FDA CFR: Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.*

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