Twin-Cath® Two-Lumen Intravenous Catheter

Rx only

Intended Purpose:

The Arrow Twin-Cath® is intended to provide short-term (<30 days) peripheral venous access.

Indications for Use:

The Arrow Twin-Cath permits peripheral venous access.

Patient Target Group:

Intended to be used in patients with anatomy suitable for use with the device and who require peripherally compatible medications.

Contraindications:

None known.

Clinical Benefits to be Expected:

The TwinCath provides short term venous access and permits administration of two medications simultaneously with just one venous puncture.

Performance Characteristics:

Allows venous access through two independent fluid pathways.

Contains Hazardous Substance:

Components manufactured using Stainless Steel can contain > 0.1% weight by weight of Cobalt (CAS # 7440-48-4) which is considered a category 1B CMR (Carcinogenic, mutagenic or toxic to reproduction) substance. The amount of Cobalt in the Stainless Steel components has been evaluated and considering the intended purpose and toxicological profile of the devices there is no biological safety risk to patients when using the devices as instructed within this IFU.

⚠️ General Warnings and Precautions

Warnings:

1. Sterile, Single use: Do not reuse, reprocess or resterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death. Reprocessing of medical devices intended for single use only may result in degraded performance or a loss of functionality.
2. Read all package insert warnings, precautions and instructions prior to use. Failure to do so may result in severe patient injury or death.
3. Air embolism can occur if air is allowed to enter a vascular access device or vein. Do not leave open needles or uncapped, unclamped catheters in venous puncture site. Use only securely tightened Luer-Lock connections with any vascular access device to guard against inadvertent disconnection.
4. Do not use scissors to remove dressing to reduce the risk of damage to the catheter.
5. Clinicians must be aware of complications/undesirable side-effects associated with peripheral intravenous catheters including, but not limited to:

- infiltration
- migration
- phlebitis
- air embolism
- catheter embolism
- bacteremia
- septicemia
- thrombosis
- inadvertent arterial puncture
- nerve injury/damage
- hematoma
- hemorrhage

Precautions:

1. Do not alter the catheter or any other kit/set component during insertion, use or removal.
2. Procedure must be performed by trained personnel well versed in anatomical landmarks, safe technique and potential complications.
3. Use standard precautions and follow institutional policies for all procedures including safe disposal of devices.
4. If the package is damaged or unintentionally opened before use do not use the device. Dispose of the device.
5. Storage conditions for these devices require that they are kept dry and out of direct sunlight.
6. Limit placement of catheter to large venous channels of upper arm and forearm.
7. Do not secure, staple and/or suture directly to outside diameter of catheter body or extension lines to reduce risk of cutting or damaging the catheter or impeding catheter flow. Secure only at indicated stabilization locations.
8. Catheter should be routinely inspected for desired flow rate, security of dressing, correct position, and proper Luer-Lock connection.
9. Maintain insertion site with regular meticulous redressing using aseptic technique.
10. For blood sampling temporarily shut off distal port if solution is being infused.
11. Some disinfectants used at catheter insertion site contain solvents which can weaken the catheter material. Alcohol, acetone, and polyethylene glycol can weaken the structure of polyurethane materials. These agents may also weaken the adhesive bond between catheter stabilization device and skin.
 - Do not use acetone on catheter surface.
 - Do not use alcohol to soak catheter surface or allow alcohol to dwell in a catheter lumen to restore catheter patency or as an infection prevention measure.
 - Do not use polyethylene glycol containing ointments at insertion site.
 - Take care when infusing drugs with a high concentration of alcohol.
 - Allow insertion site to dry completely prior to skin puncture and before applying dressing.
 - Do not allow kit components to come into contact with alcohol.

Kits/Sets may not contain all accessory components detailed in these instructions for use. Become familiar with instructions for individual component(s) before beginning the procedure.

A Suggested Procedure: Use Sterile Technique.

1. Prepare puncture site as required.
2. Prepare catheter for insertion. Flush proximal port through injection cap. Leave guard in place, hold catheter upright, and flush. This will activate the catheter's hydrophilic coating.

⚠️ Precaution: Do not allow flush solution to go beyond tip of catheter to reduce interference with blood flashback.

3. Puncture vessel using a continuous, controlled, slow forward motion. Avoid transfixing both vessel walls. Blood flashback in hub of introducer needle indicates successful entry into vessel. Aspiration may be required.
4. Advance catheter and needle, as a unit, approximately one centimeter, to ensure that vessel dilation is complete.
5. Hold needle in position and advance catheter forward into vessel. Remove introducer needle.

⚠️ Precaution: Do not reinsert needle into catheter to reduce the risk of catheter damage.

6. Attach stopcock, injection cap, or connecting tubing to distal hub.

⚠️ Precaution: To lessen the risk of disconnects, only securely tightened Luer-Lock connections should be used with this device.

⚠️ Precaution: Do not begin infusion until proximal lumen placement is verified.

7. Verify proximal lumen placement by aspirating blood from the proximal port through the extension line and flush.
8. Attach proximal hub to connecting line or heparinlock through injection cap using standard hospital protocol. A clamp is provided to occlude flow through the proximal lumen.

⚠️ Precaution: Open clamp prior to infusion.

9. Secure catheter to patient. Suture wing or suture ring on juncture hub of catheter may be used.

⚠️ Precaution: If using suture ring, suture should be snug but not tight. Tightening will constrict lumens.

10. Dress puncture site per hospital protocol.

⚠️ Precaution: Catheter should be routinely inspected for desired flow rate, security of dressing, correct position, and proper Luer-Lock connection.

Infusate Guidelines

1. Two diluted medications, continuous or intermittent, can be administered simultaneously.
2. One undiluted medication and one diluted medication can be administered simultaneously.

⚠️ Warning: Simultaneous administration of two undiluted medications is not recommended since there is clear *in-vitro* evidence of medication interaction at normal peripheral venous flow rates.

Care and Maintenance:

Dressing:

Dress according to institutional policies, procedures, and practice guidelines. Change immediately if the integrity becomes compromised e.g. dressing becomes damp, soiled, loosened or no longer occlusive.

Catheter Patency:

Maintain catheter patency according to institutional policies, procedures and practice guidelines.

Catheter Removal Instructions:

1. Remove dressing.
2. Remove from catheter securement device.
3. Remove catheter by slowly pulling it parallel to skin. If resistance is met while removing catheter **STOP**

⚠️ Precaution: Catheter should not be forcibly removed, doing so may result in catheter breakage and embolization. Follow institutional policies and procedures for difficult to remove catheter.

4. Apply direct pressure to site until hemostasis is achieved followed by an ointment-based occlusive dressing.

⚠️ Warning: Residual catheter track remains an air entry point until site is epithelialized. Occlusive dressing should remain in place for at least 24 hours or until site appears epithelialized.

5. Document catheter removal procedure including confirmation that entire catheter length has been removed per institutional policies and procedures.

For reference literature concerning patient assessment, clinician education, insertion technique, and potential complications associated with this procedure, consult standard textbooks, medical literature, and Arrow International LLC website: www.teleflex.com

A pdf copy of this IFU is located at www.teleflex.com/IFU

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority. The contacts of national competent authorities (Vigilance Contact Points) and further information can be found on the following European Commission website: https://ec.europa.eu/growth/sectors/medical-devices/contacts_en

en

Symbol Glossary: Symbols are in compliance with ISO 15223-1. Some symbols may not apply to this product. Refer to product labeling for symbols that apply specifically to this product.

								
Caution	Medical device	Consult instructions for use	Contains hazardous substances	Do not reuse	Do not resterilize	Sterilized by ethylene oxide	Single sterile barrier system	Keep away from sunlight
								
Keep dry	Do not use if package is damaged	Not made with natural rubber latex	Catalogue number	Lot number	Use by	Manufacturer	Date of manufacture	Importer

Arrow, the Arrow logo, Teleflex, the Teleflex logo and Twin-Cath are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. © 2023 Teleflex Incorporated. All rights reserved.

"Rx only" is used within this labeling to communicate the following statement as presented in the FDA CFR : Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.

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