

REF **CDC-21242-1A**

9 Fr. distal lumen | **2** Lumen | **11.5** cm device length | **.035** inch dia. spring-wire guide

ErgoPack®

- 1: Two-Lumen Central Venous Access Device¹: 115 mm Radiopaque Polyurethane with 9 Fr. Distal Lumen, ARROWgard[®] Antimicrobial Surface Treatment², Integral Hemostasis Valve, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Injection Needle: SafetyGlide™ 25 Ga. x 1" (2.54 cm)
- 1: Local Anesthetic Component
- 1: Pressure Transduction Probe
- 1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow[®] Raulerson Spring-Wire Introduction Syringe
- 1: Injection Needle: SafetyGlide™ 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
- 1: Syringe: 3 mL Luer-Lock

- 1: Tissue Dilator
- 1: 5 mL 1% Lidocaine HCl Solution
- 1: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®4 One-Step Solution with Hi-Lite Orange™ Tint (non-sterile solution)
- 1: Obturator: 8 Fr.
- 1: Cath-Gard® Catheter Contamination Shield: 80 cm with TwistLock™ Adapter
- 1: Dust Cap: Vented
- 1: Dust Cap: Non-Vented
- 1: SharpsAway® Disposal Cup
- 1: SharpsAway II™ Locking Disposal Cup
- 1: Maximal Barrier Drape™ with 4" fenestration
- 1: Needle Holder
- 1: Safety Scalpel: #11
- 1: Stopcock: 4-way High-Flow

- 1: Checklist/CLIP Sheet
- 1: Sterile Procedure Sign
- 1: Medication Label: 1% Lidocaine
- 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 5: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Surgical Apparel: Impervious Gown
- 1: Surgical Apparel: Bouffant Cap
- 1: Surgical Apparel: Mask with Eye Shield
- 1: Suture: 3-0 Silk with Curved Needle
- 1: Tubing*: 20" (50 cm) Extension
- 1: HemoHopper® Fluid Receptacle

¹Licensed under US Patent No. 5,250,038.²Licensed under US Patent No. 6,706,024.


³A trademark of Becton, Dickinson and Company.

⁴A registered trademark of CareFusion or one of its subsidiaries.

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Contraindications: The ARROW[®]g[®]ard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine acetate, silver sulfadiazine and/or sulfa drugs.

Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.



* Priming volumes are approximate and are done without accessories.
† Based on 100 cm head height.

