



## ARROWgard Blue PLUS<sup>®</sup> Pressure Injectable Two-Lumen CVC Kit

REF CDC-42802-P1A

8 Fr. | 2 Lumen | 16 cm catheter length | .032 inch dia. spring-wire guide

**Contents:**

- 1: Two-Lumen Indwelling Catheter: 8 Fr. x 16 cm Pressure Injectable Radiopaque Polyurethane with Blue FlexTip<sup>®</sup>, ARROWgard Blue PLUS<sup>®</sup> Antimicrobial Surface Treatment<sup>1</sup>, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Injection Needle: SafetyGlide<sup>™</sup> 25 Ga. x 1" (2.54 cm)
- 1: Local Anesthetic Component
- 1: Pressure Transduction Probe
- 1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow<sup>®</sup> Raulerson Spring-Wire Introduction Syringe
- 1: Injection Needle: SafetyGlide<sup>™</sup> 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
- 1: Syringe: 3 mL Luer-Lock
- 1: Tissue Dilator: 9 Fr. (3.0 mm) x 10.2 cm
- 1: 5 mL 1% Lidocaine HCl Solution
- 1: 3 mL Applicator 2% CHG and 70% IPA Chloraprep<sup>®</sup> One-Step Solution with Hi-Lite Orange<sup>™</sup> Tint (non-sterile solution)
- 1: SecondSite<sup>™</sup> Adjustable Hub: Fastener
- 1: SecondSite<sup>™</sup> Adjustable Hub: Catheter Clamp
- 1: SharpsAway II<sup>™</sup> Locking Disposal Cup

- 1: SharpsAway<sup>®</sup> Disposal Cup
- 1: Maximal Barrier Drape<sup>™</sup> with 4" fenestration
- 1: Needle Holder
- 1: Safety Scalpel: #11
- 1: Flow Rate Information Card with Injection Log
- 1: Sterile Procedure Sign
- 1: Medication Label: 1% Lidocaine
- 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 5: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Surgical Apparel: Impervious Gown
- 1: Surgical Apparel: Bouffant Cap
- 1: Surgical Apparel: Mask with Eye Shield
- 1: Suture: 3-0 Silk with Curved Needle

<sup>1</sup>Licensed under US Patent Nos. 6,706,024; 6,872,195 and 7,329,412.  
<sup>2</sup>A trademark of Becton, Dickinson and Company.  
<sup>3</sup>A registered trademark of CareFusion or one of its subsidiaries.



STERILE EO



Not made with natural rubber latex



Store between 20 - 25°C (68 - 77°F).

Rx only. **Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.**

Sterilized using ethylene oxide unless otherwise indicated in the contents list.


**California Prop. 65**

**Warning: Cancer and Reproductive Harm [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).**

Fluid path components are non-pyrogenic.

**Contraindications: The ARROWgard Blue Plus antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).**

Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)	MAX Pressure Injection Flow Rate‡ (mL/sec)
Distal (14 Ga.)	0.70	7650	10
Proximal (14 Ga.)	0.71	6000	10

\* Priming volumes are approximate and are done without injection caps.  
 † Based on 100 cm head height.  
 ‡ Pressure injection flow rates are determined using media of 11.8 centipoise viscosity, with 60" administration set tubing.

Manufacturer:  
 Arrow International, Inc.  
 Subsidiary of Teleflex Incorporated  
 2400 Bernville Road  
 Reading, PA 19605 USA

**Teleflex<sup>®</sup>**  
[www.teleflex.com](http://www.teleflex.com)

C-42802-101F (04/18)



## ARROWgard Blue PLUS<sup>®</sup> Pressure Injectable Two-Lumen CVC Kit

REF CDC-42802-P1A

8 Fr. | 2 Lumen | 16 cm catheter length | .032 inch dia. spring-wire guide