







## **ARROW ENDURANCE™ Extended Dwell Peripheral Catheter System**



8 cm catheter length

inch dia.

- Catheter System: 18 Ga. (1.30 mm) x 80 mm (3.15") ARROW ENDURANCE™ over 20 Ga. (1.0 mm) HW Introducer Needle with Integral 0.010" (0.26 mm) Straight Soft Tip Wire Guide
- Skin Protectant Prep
- 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®¹ One-Step Solution with Hi-Lite Orange™ Tint Towel: 24" x 36" (61 cm x 91 cm)
  Arm Drape: 36" x 46" with 3" x 5" Fenestration

- Label: Venous/Arterial Placement Paper Tape Measure
- Tourniquet

- Tourniquet
   Gauze Pad: 2" x 2" (5 cm x 5 cm)
   Gauze Pad: 4" x 4" (10 cm x 10 cm)
   Dressing: Tegaderm<sup>®</sup>2 10 cm x 12 cm
   Dressing: STATLOCK<sup>®</sup>3 Catheter Stabilization Device
   Surgical Apparel: Bouffant Cap
   Surgical Apparel: Mask

 $\mathsf{Rx}$  only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

⚠ California Prop. 65
Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.

Fluid path components are non-pyrogenic.

Priming Volume	Gravity Flow Rate	Pump Flow Rate	MAX Injector Pressure Limit	MAX Pressure Injection Flow Rate per Injectate Viscosity		
volunic				4.7 cP	11.8 cP	26.6 cP
0.42 mL	3682 mL/hr	11782 mL/hr	325 psi	10.0 mL/sec	8.0 mL/sec	4.5 mL/sec

Consult instructions for use for details on the testing performed to generate these values.





















Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA



EDC-00818-201B (2021-04)







ARROW ENDURANCETM Extended Dwell Peripheral Catheter System  $18^{\text{Ga.}}$   $8^{\text{cm}}_{\text{length}}$  .010 wire guide



<sup>&</sup>lt;sup>1</sup>A registered trademark of CareFusion or one of its subsidiaries. <sup>2</sup>A registered trademark of 3M Company. <sup>3</sup>A registered trademark of C. R. Bard, Inc.