



Arterial Catheterization Kit

REF **NA-04550-1A**

20 Ga. | **12** cm catheter length | **.025** inch dia. spring-wire guide

Contents:

- | | | |
|--|---|---|
| 1: Indwelling Catheter: 20 Ga. x 12 cm Radiopaque Polyurethane with Integral Suture Wing | 2: 3 mL Applicator 2% CHG and 70% IPA Chloraprep® ² One-Step Solution with Hi-Lite Orange™ Tint (non-sterile solution) | 1: Surgical Apparel: Bouffant Cap |
| 1: Spring-Wire Guide: .025" (0.64 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer | 1: SharpsAway® Disposal Cup | 1: Surgical Apparel: Mask with Eye Shield |
| 1: Catheter: 20 Ga. x 1-3/4" (4.45 cm) Radiopaque over 22 Ga. TW Introducer Needle | 1: SharpsAway II™ Locking Disposal Cup | 1: Suture: 3-0 Silk with Curved Needle |
| 1: Introducer Needle: 20 Ga. x 2" (5.08 cm) XTW | 1: Drape: 41" x 55" with 4" fenestration | ¹ A trademark of Becton, Dickinson and Company. |
| 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) RW | 1: Towel: 13" x 18" (33 cm x 46 cm) | ² A registered trademark of CareFusion or one of its subsidiaries. |
| 1: Injection Needle: SafetyGlide™ ¹ 25 Ga. x 1" (2.54 cm) | 1: Needle Holder | ³ A registered trademark of 3M Company. |
| 1: Injection Needle: SafetyGlide™ ¹ 22 Ga. x 1-1/2" (3.81 cm) RW | 1: Safety Scalpel: #11 | |
| 1: Local Anesthetic Component | 1: Medication Label: 1% Lidocaine | |
| 1: Syringe: 3 mL Luer-Lock | 4: Gauze Pad: 2" x 2" (5 cm x 5 cm) | |
| 1: 5 mL 1% Lidocaine HCl Solution | 4: Gauze Pad: 4" x 4" (10 cm x 10 cm) | |
| | 1: Dressing: Tegaderm® ³ 10 cm x 12 cm | |
| | 1: Surgical Apparel: Impervious Gown | |

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.



Not made with natural rubber latex



Store between 20-25°C (68-77°F).