





## **Arterial Line Kit**

20<sup>Ga.</sup> 12<sup>cm</sup> catheter length 0

- Indwelling Catheter: 20 Ga. x 12 cm Radiopaque FEP Spring-Wire Guide: .025" (0.64 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on One End "J" Tip on Other) Injection Needle: 25 Ga. x 1" (2.54 cm) Injection Needle: 18 Ga. x 1-1/2" (3.81 cm) Introducer Needle: 20 Ga. x 2" (5.08 cm) XTW Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm) Syringe: 3 ml Just Jock

- Injection Needle: 18 Ga. x 1-1/2" (3.81 cm)
   Introducer Needle: 20 Ga. x 2" (5.08 cm) XTW
   Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
   Syringe: 3 mL Luer-Lock
   5 mL 1% Lidocaine HCI Solution and Alcohol Prep
   3 mL Applicator 2% CHG and 70% IPA ChloraPrep®¹ One-Step Solution with Hi-Lite Orange™ Tint
   Drape: 18" x 26" (46 cm x 66 cm) with 3" (7.6 cm) fenestration
   Scalpel: #11
   Gauze Pad: 2" x 2" (5 cm x 5 cm)
   Suture: 3-0 Braided Silk with Straight Needle

<sup>1</sup>A registered trademark of CareFusion or one of its subsidiaries.

Rx only. Warning: Read all package insert warnings, cautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

⚠ California Prop. 65
Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.

Fluid path components are non-pyrogenic.





















Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA



IK-04510-105D (2022-06)







