





Arterial Catheterization Kit

20^{Ga.} 12^{cm} catheter length 0

- Indwelling Catheter: 20 Ga. x 12 cm Radiopaque FEP Spring-Wire Guide: .025" (0.64 mm) dia. x 13-1/8" (33 cm) (Straight Soft Tip on Both Ends) Catheter: 20 Ga. x 1-1/4" (3.2 cm) Jelco®1 Catheter over 22 Ga. Introducer Needle

- 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)
 1: Introducer Needle: 20 Ga. x 2" (5.08 cm) XTW
 2: Syringe: 3 mL Luer-Lock
 1: Drape: 18" x 26" (46 cm x 66 cm) with 3" (7.6 cm) fenestration
- 1: Prep Sponge Swab
 3: Gauze Pad: 2" x 2" (5 cm x 5 cm)

¹A registered trademark of Medex, Inc.

Rx only. Warning: Read all package insert warnings, cautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

\(\triangle \) California Prop. 65

Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.

Fluid path components are non-pyrogenic.



















Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA



IS-04100-105A (2020/04)







