





7 Fr. **3** Lumen

20 cm catheter length .025 inch dia. spring-wire guide

Arrowg+ard Blue Plus® Three-Lumen CVC

- Arrowg+ard Blue Plus®¹ Three-Lumen Catheter: 7 Fr. (2.5 mm OD) x 20 cm Spring-Wire Guide, Marked: .025" (0.64 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End "J" Tip on Other) with Arrow Advancer
- Catheter: 20 Ga. x 1-3/4" (4.45 cm) Radiopaque over 22 Ga. TW Introducer
- Needle: 18 Ga. x 2-1/2" (6.35 cm) RW Injection Needle: SafetyGlide™2 25 Ga. x 5/8" (1.60 cm) RW Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
- Pressure Transduction Probe
- Introducer Needle: 20 Ga. x 1-1/2" (3.81 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe Injection Needle: SafetyGlide™2 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip
- Syringe Syringe: 3 mL Luer-Lock
- Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
- Issue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
 5 mL 19k Lidocaine HCl Solution and Alcohol Prep
 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®3 One-Step Solution with Hi-Lite Orange™ Tint
 SecondSite™ Adjustable Hub: Fastener
 SecondSite™ Adjustable Hub: Catheter Clamp
 SharpsAway® II Locking Disposal Cup

- SharpsAway® Disposal Cup Drape: 41" x 55" with 4" fenestration

- Drape: 41" x 55" with 4" fenestration
 Safety Scalpel: #11
 Gauze Pad: 2" x 2" (5 cm x 5 cm)
 Gauze Pad: 4" x 4" (10 cm x 10 cm)
 Dressing: STATLOCK®4 CV Plus Catheter Stabilization Device
 Suture: 3-0 Braided Silk with Straight Needle
- ¹Licensed under US Patent No. 7.329.412.
- ³A registered trademark of CareFusion or one of its subsidiaries. ⁴A registered trademark of C. R. Bard, Inc.

Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

California Prop. 65

WARNING: Cancer and Reproductive Harm.
 www.P65Warnings.ca.gov
 Not made with natural rubber latex.

Store between 20 - 25°C (68 - 77°F)

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.

Contraindications: The Arrowg+ard Blue Plus antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs.

Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)	
Distal (16 Ga.)	0.39	2283	
Medial (18 Ga.)	0.39	913	
Proximal (18 Ga.)	0.40	996	

^{*} Priming volumes are approximate and are done without accessories.









Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA



LBL051352 R00 (2021-11)











7 Fr. 3 Lumen 20 cm catheter length length guide



[†] Flow rate values are approximate and are determined using deionized water at 100 cm head height.