



REF CDC-25142-X1A

14 Fr. 2 Lumen 20 cm catheter length .035 inch dia. spring-wire guide

ARROW

ErgoPack®

# Arrowg+ard Blue® Two-Lumen Hemodialysis Catheterization for High Volume Infusions

## Contents:

- |  |   |
|--|---|
| <ul style="list-style-type: none"> <li>1: Two-Lumen Catheter: 14 Fr. (4.72 mm OD) x 20 cm Radiopaque Polyurethane with Blue FlexTip®, Arrowg+ard® Antimicrobial Surface Treatment, Extension Line Clamps</li> <li>1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer</li> <li>1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle</li> <li>1: Injection Needle: SafetyGlide™<sup>1</sup> 25 Ga. x 1" (2.54 cm)</li> <li>1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)</li> <li>1: Pressure Transduction Probe</li> <li>1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe</li> <li>1: Injection Needle: SafetyGlide™<sup>1</sup> 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe</li> <li>1: Syringe: 3 mL Luer-Lock</li> <li>1: Tissue Dilator: 14 Fr. (4.7 mm) x 14 cm</li> <li>1: Tissue Dilator: 12 Fr. (4.0 mm) x 14 cm</li> <li>2: 10 mL Luer-Lock Pre-Filled Saline Syringe</li> <li>1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep</li> <li>1: 3 mL Applicator 2% CHG and 70% IPA Chloraprep®<sup>2</sup> One-Step Solution with Hi-Lite Orange™ Tint</li> </ul> | <ul style="list-style-type: none"> <li>1: SharpsAway® Disposal Cup</li> <li>1: SharpsAway® II Locking Disposal Cup</li> <li>1: Maximal Barrier Drape™ with 4" fenestration</li> <li>1: Towel</li> <li>1: Needle Holder</li> <li>1: Safety Scalpel: #11</li> <li>1: Checklist/CLIP Sheet</li> <li>1: Sterile Procedure Sign</li> <li>1: Medication Label: 1% Lidocaine</li> <li>2: Gauze Pad: 2" x 2" (5 cm x 5 cm)</li> <li>5: Gauze Pad: 4" x 4" (10 cm x 10 cm)</li> <li>1: Dressing: Tegaderm®<sup>3</sup> 10 cm x 12 cm</li> <li>1: Surgical Apparel: Impervious Gown</li> <li>1: Dressing: BIOPATCH®<sup>4</sup></li> <li>1: Surgical Apparel: Bouffant Cap</li> <li>1: Surgical Apparel: Mask with Eye Shield</li> <li>1: Suture: 3-0 Braided Silk with Straight Needle</li> <li>1: Suture: 3-0 Silk with Curved Needle</li> <li>1: HemoHopper® Fluid Receptacle</li> </ul> |
|--|---|

<sup>3</sup>A registered trademark of 3M Company.  
<sup>4</sup>A registered trademark of Johnson & Johnson Corporation.

## Not made with natural rubber latex.


Store between 20 – 25°C (68 – 77°F).  
 Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.

**Contraindications:** The Arrow large-bore two-lumen catheter is not designed for long-term (≥ 30 days) hemodialysis or for use in patients with thrombosed vessels.  
 The Arrowg+ard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine, and/or sulfa drugs.



Priming Volume*	
Distal	1.70 mL
Proximal	1.50 mL

\*Priming volumes are approximate and are done without accessories.



† In vitro tests performed using 3:1 glycerine/saline at 37°C.

