



REF SB-09803

8.5 Fr. 10 cm sheath length .025 inch dia. spring-wire guide

PSI Kit for use with 7 - 7.5 Fr. Catheters

Contents:

- | | |
|--|--|
| 1: Arrow-Flex® Sheath: 8.5 Fr. (2.8 mm tip ID) x 10 cm and Tissue Dilator | 2: Prep Sponge Swab |
| 1: Spring-Wire Guide, Marked: .025" (0.64 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer | 2: Gauze Pad: 2" x 2" (5 cm x 5 cm) |
| 1: Catheter: 20 Ga. x 1-3/4" (4.45 cm) Radiopaque over 22 Ga. TW Introducer Needle | 5: Gauze Pad: 4" x 4" (10 cm x 10 cm) |
| 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) | 1: Suture: 3-0 Braided Silk with Straight Needle |
| 1: Injection Needle: 25 Ga. x 1" (2.54 cm) | 1: Tray: Prep |
| 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) RW | |
| 1: Pressure Transduction Probe | |
| 1: Introducer Needle: 20 Ga. x 1-1/2" (3.81 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe | |
| 1: Syringe: 3 mL Luer-Lock | |
| 1: Syringe: 5 mL Luer-Slip | |
| 1: Obturator: 8 Fr. | |
| 1: Cath-Gard® Catheter Contamination Shield: 80 cm with Sterile Tape | |
| 1: SharpsAway® Disposal Cup | |
| 1: Drape: 24" x 36" (61 cm x 91 cm) with 4" (10 cm) fenestration | |
| 1: Drape: 24" x 36" (61 cm x 91 cm) | |
| 1: Scalpel: #11 | |
| 1: Stopcock: 3-way High-Flow | |

Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Caution: This Product Contains Dry Natural Rubber.

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.



LOT



Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

LBL058341 R01 (2022-02)



REF SB-09803

8.5 Fr. 10 cm sheath length .025 inch dia. spring-wire guide

PSI Kit for use with 7 - 7.5 Fr. Catheters