



REF CDC-29903-X1A

(IPNXXXXXX)

9^{Fr.}10^{cm}
sheath
length.035^{inch dia.}
spring-wire
guide

ARROW

ErgoPack®

Arrowgard Blue® PSI for use with 7.5 - 8 Fr. Catheters

**Contents:**

- | | |
|--|--|
| 1: Arrowgard Blue® Arrow-Flex® Sheath: 9 Fr. (2.9 mm tip ID) x 10 cm and Tissue Dilator | 1: SharpsAway® II Locking Disposal Cup |
| 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer | 1: SharpsAway® Disposal Cup |
| 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle | 1: Maximal Barrier Drape™ with 4" fenestration |
| 1: Injection Needle: SafetyGlide™ ¹ 25 Ga. x 1" (2.54 cm) | 1: Towel |
| 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm) | 1: Needle Holder |
| 1: Pressure Transduction Probe | 1: Safety Scalpel: #11 |
| 1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe | 1: Stopcock: 4-way High-Flow |
| 1: Injection Needle: SafetyGlide™ ¹ 23 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe | 1: Checklist/CLIP Sheet |
| 1: Syringe: 3 mL Luer-Lock | 1: Sterile Procedure Sign |
| 2: 10 mL Luer-Lock Pre-Filled Saline Syringe | 1: Medication Label: 1% Lidocaine |
| 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep | 2: Gauze Pad: 2" x 2" (5 cm x 5 cm) |
| 1: 3 mL Applicator 2% CHG and 70% IPA Chloraprep® ² One-Step Solution with HI-Lite Orange™ Tint | 5: Gauze Pad: 4" x 4" (10 cm x 10 cm) |
| 1: Obturator: 8 Fr. | 1: Dressing: Tegaderm® ³ 10 cm x 12 cm |
| 1: Cath-Gard® Catheter Contamination Shield: 80 cm with TwistLock™ Adapter | 1: Surgical Apparel: Impervious Gown |
| 1: Dust Cap: Non-Vented | 1: Dressing: BIOPATCH® ⁴ |
| | 1: Surgical Apparel: Bouffant Cap |
| | 1: Surgical Apparel: Mask with Eye Shield |
| | 1: Suture: 3-0 Silk with Curved Needle |
| | 1: Tubing±: 20" (50 cm) Extension |
| | 1: HemoHopper® Fluid Receptacle |
| | ¹ A trademark of Becton, Dickinson and Company. |

²A registered trademark of CareFusion or one of its subsidiaries.³A registered trademark of 3M Company.⁴A registered trademark of Johnson & Johnson Corporation.**±Precaution: Contains Phthalate: DEHP.**

Warning: The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction. Proceeding from the present state of scientific knowledge, risks for male premature infants cannot be excluded in the case of long-term exposure or application. Medical products containing phthalates should be used only temporarily with pregnant women, nursing mothers, babies and infants.

Not made with natural rubber latex.

Store between 20 – 25°C (68 – 77°F).

Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.

Contraindications: The Arrowgard Blue antimicrobial sheath introducer is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine, and/or sulfa drugs (refer to product instructions for references).

