



REF CDA-29803-1A

(IPNXXXXXX)

8.5 Fr.

10 cm
sheath
length.035 inch dia.
spring-wire
guide

ARROW

ErgoPack®

Arrowgard Blue® PSI for use with 7.5 - 8 Fr. Catheters

**Contents:**

- | | |
|--|---|
| 1: Arrowgard Blue® Arrow-Flex® Sheath: 8.5 Fr. (2.8 mm tip ID) x 10 cm and Tissue Dilator | 1: Forceps |
| 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer | 1: Safety Scalpel: #11 |
| 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle | 1: Stopcock: 4-way High-Flow |
| 1: Pressure Transduction Probe | 1: Checklist/CLIP Sheet |
| 1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe | 1: Sterile Procedure Sign |
| 1: Injection Needle: Eclipse™ ¹ 23 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe | 2: Gauze Pad: 2" x 2" (5 cm x 5 cm) |
| 1: Injection Needle: Eclipse™ ¹ 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock Syringe | 5: Gauze Pad: 4" x 4" (10 cm x 10 cm) |
| 1: Syringe: 10 mL Luer-Lock | 1: Dressing: Tegaderm® ² 10 cm x 12 cm |
| 1: Obturator: 8 Fr. | 1: Surgical Apparel: Impervious Gown |
| 1: Cath-Gard® Catheter Contamination Shield: 80 cm with TwistLock™ Adapter | 1: Surgical Apparel: Bouffant Cap |
| 1: Dust Cap: Non-Vented | 1: Surgical Apparel: Mask with Eye Shield |
| 1: SharpsAway® II Locking Disposal Cup | 1: Suture: 3-0 Silk with Curved Needle |
| 1: SharpsAway® Disposal Cup | 1: Tubing: 20" (50 cm) Extension |
| 1: Maximal Barrier Drape™ with 4" fenestration | 1: HemoHopper® Fluid Receptacle |
| 1: Towel | |
| 1: Needle Holder | |

¹A trademark of Becton, Dickinson and Company.²A registered trademark of 3M Company.**±Precaution: Contains Phthalate: DEHP.**

Warning: The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction. Proceeding from the present state of scientific knowledge, risks for male premature infants cannot be excluded in the case of long-term exposure or application. Medical products containing phthalates should be used only temporarily with pregnant women, nursing mothers, babies and infants.

Not made with natural rubber latex.

Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F).

Contraindications: The Arrowgard Blue antimicrobial sheath introducer is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine, and/or sulfa drugs (refer to product instructions for references).

