



REF CDC-29803-XN1A

(IPNXXXXXX)

8.5 Fr.

10 cm  
sheath  
length.035 inch dia.  
spring-wire  
guide

ARROW

ErgoPack.  
Complete

# Arrowgard Blue<sup>®</sup> PSI for use with 7.5 - 8 Fr. Catheters

UDI

**Contents:**

- |   |   |
|---|---|
| 1: Arrowgard Blue <sup>®</sup> Arrow-Flex <sup>®</sup> Sheath: 8.5 Fr. (2.8 mm tip ID) x 10 cm and Tissue Dilator   | 1: Towel                                    |
| 1: Spring-Wire Guide, Nitinol, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow <sup>®</sup> GlideWheel <sup>™</sup> Wire Advancer | 1: Needle Holder                            |
| 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle  | 1: Safety Scalpel: #11                      |
| 1: Injection Needle: SafetyGlide <sup>™</sup> 1 25 Ga. x 1" (2.54 cm)   | 1: Stopcock: 4-way High-Flow                |
| 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)   | 1: Checklist/CLIP Sheet                     |
| 1: Pressure Transduction Probe  | 1: Sterile Procedure Sign                   |
| 1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow <sup>®</sup> Raulerson Spring-Wire Introduction Syringe  | 1: Medication Label: 1% Lidocaine           |
| 1: Injection Needle: SafetyGlide <sup>™</sup> 1 23 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe  | 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)         |
| 1: Syringe: 3 mL Luer-Lock  | 5: Gauze Pad: 4" x 4" (10 cm x 10 cm)       |
| 2: 10 mL Luer-Lock Pre-Filled Saline Syringe  | 1: Surgical Apparel: Impervious Gown        |
| 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep  | 1: Dressing: I.V. Advanced 8.5 cm x 11.5 cm |
| 2: 3 mL Applicator 2% CHG and 70% IPA Chloraprep <sup>®2</sup> One-Step Solution with Hi-Lite Orange <sup>™</sup> Tint  | 1: Dressing: BIOPATCH <sup>®3</sup>         |
| 1: Obturator: 8 Fr.   | 1: Transducer Cover: 14 cm x 147 cm         |
| 1: Cath-Gard <sup>®</sup> Catheter Contamination Shield: 80 cm with TwistLock <sup>™</sup> Adapter  | 1: Surgical Apparel: Mask with Eye Shield   |
| 1: Dust Cap: Non-Vented   | 1: Surgical Apparel: Bouffant Cap           |
| 1: SharpsAway <sup>®</sup> II Locking Disposal Cup  | 1: Suture: 3-0 Silk with Curved Needle      |
| 1: SharpsAway <sup>®</sup> Disposal Cup   | 1: Tubing±: 20" (50 cm) Extension           |
| 1: Maximal Barrier Drape <sup>™</sup> with 4" fenestration  | 1: HemoHopper <sup>®</sup> Fluid Receptacle |

<sup>1</sup>A trademark of Becton, Dickinson and Company.<sup>2</sup>A registered trademark of CareFusion or one of its subsidiaries.<sup>3</sup>A registered trademark of Johnson & Johnson Corporation.**±Precaution: Contains Phthalate: DEHP.**

**Warning: The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction. Proceeding from the present state of scientific knowledge, risks for male premature infants cannot be excluded in the case of long-term exposure or application. Medical products containing phthalates should be used only temporarily with pregnant women, nursing mothers, babies and infants.**

**Not made with natural rubber latex.**

Store between 20 – 25°C (68 – 77°F).

Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.

**Contraindications: The Arrowgard Blue antimicrobial sheath introducer is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine, and/or sulfa drugs (refer to product instructions for references).**

