



REF ASK-04020-MHS

20 Ga.

4.45 cm catheter length

.018 inch dia. spring-wire guide

Radial Artery Catheterization

Contents:

- 1: Catheter: 20 Ga. x 4.45 cm Radiopaque Polyurethane over 22 Ga. TW
Introducer Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide

1: Injection Needle: 18 Ga. x 1-1/2" (3.81 cm)

1: Injection Needle: 27 Ga. x 1/2" (1.27 cm)

1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)

1: Syringe: 3 mL Luer-Lock

1: Skin Protectant Prep

1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep

1: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®¹ One-Step Solution with Hi-Lite Orange™ Tint

1: Wing Clip: Suture

1: SharpsAway® II Locking Disposal Cup

1: SharpsAway® Disposal Cup

1: Towel: 13" x 18" (33 cm x 46 cm)

1: Drape: 18" x 26" with 2" x 3" fenestration, with adhesive

1: Safety Scalpel: #11

1: Bendable Armboard‡

4: Gauze Pad: 2" x 2" (5 cm x 5 cm)

4: Gauze Pad: 4" x 4" (10 cm x 10 cm)

1: Dressing: Tegaderm®² 10 cm x 12 cm
- 1: Suture: 3-0 Braided Silk with Straight Needle

¹A registered trademark of CareFusion or one of its subsidiaries.

²A registered trademark of 3M Company.

Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU

California Prop. 65

⚠ **WARNING:** Cancer and Reproductive Harm.
www.P65Warnings.ca.gov

‡MR Conditional

Not made with natural rubber latex.

Store between 20 – 25°C (68 – 77°F).

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.



Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

LBL059310 R03 (2023-08)



REF ASK-04020-MHS

20 Ga.

4.45 cm catheter length

.018 inch dia. spring-wire guide

Radial Artery Catheterization