



REF NA-04020-1A (IPNXXXXXX)

20 Ga. 4.45 cm catheter length .018 inch dia. spring-wire guide

Arterial Catheterization

Contents:

- 1: Catheter: 20 Ga. x 4.45 cm Radiopaque Polyurethane over 22 Ga. TW Introducer Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide origin:xxxxxxxxxxxxxxxxxxxx
- 1: Injection Needle: SafetyGlide™ 25 Ga. x 1" (2.54 cm) origin:xxxxxxxxxxxxxxxxxxxx
- 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm) origin:xxxxxxxxxxxxxxxxxxxx
- 1: Syringe: 3 mL Luer-Lock origin:xxxxxxxxxxxxxxxxxxxx
- 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep origin:xxxxxxxxxxxxxxxxxxxx
- 1: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®2 One-Step Solution with Hi-Lite Orange™ Tint origin:xxxxxxxxxxxxxxxxxxxx
- 1: Wing Clip: Suture origin:xxxxxxxxxxxxxxxxxxxx
- 1: SharpsAway® Disposal Cup origin:xxxxxxxxxxxxxxxxxxxx
- 1: SharpsAway® II Locking Disposal Cup origin:xxxxxxxxxxxxxxxxxxxx
- 1: Drape: 36" x 41" (91 cm x 104 cm) with 3" x 5" (8 cm x 13 cm) fenestration origin:xxxxxxxxxxxxxxxxxxxx
- 1: Towel: 13" x 18" (33 cm x 46 cm) origin:xxxxxxxxxxxxxxxxxxxx
- 1: Needle Holder origin:xxxxxxxxxxxxxxxxxxxx
- 1: Safety Scalpel: #11 origin:xxxxxxxxxxxxxxxxxxxx
- 1: Medication Label: 1% Lidocaine origin:xxxxxxxxxxxxxxxxxxxx
- 4: Gauze Pad: 2" x 2" (5 cm x 5 cm) origin:xxxxxxxxxxxxxxxxxxxx
- 4: Gauze Pad: 4" x 4" (10 cm x 10 cm) origin:xxxxxxxxxxxxxxxxxxxx
- 1: Dressing: Tegaderm®3 10 cm x 12 cm origin:xxxxxxxxxxxxxxxxxxxx
- 1: Surgical Apparel: Bouffant Cap origin:xxxxxxxxxxxxxxxxxxxx
- 1: Surgical Apparel: Mask with Eye Shield origin:xxxxxxxxxxxxxxxxxxxx
- 1: Suture: 3-0 Silk with Curved Needle origin:xxxxxxxxxxxxxxxxxxxx

¹A trademark of Becton, Dickinson and Company.
²A registered trademark of CareFusion or one of its subsidiaries.
³A registered trademark of 3M Company.

Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU

California Prop. 65

WARNING: Cancer and Reproductive Harm. www.P65Warnings.ca.gov

Not made with natural rubber latex.

Store between 20 – 25°C (68 – 77°F).

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.

Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.



LOT LotNumber

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)LotNumber

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Set Product of xxx
Packaged in xxx



LBL059313 R01 (2024-04)

REF NA-04020-1A

20 Ga. 4.45 cm catheter length .018 inch dia. spring-wire guide



Arterial Catheterization

