



20 Ga. **3.81** cm catheter length **.018** inch dia. spring-wire guide

Contents:

- | | |
|---|---|
| 1: QuickFlash® Radial Artery Catheter: 20 Ga. x 1-1/2" (3.81 cm) Radiopaque Polyurethane over 21 Ga. TW Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide | 1: Surgical Apparel: Mask with Eye Shield |
| 1: Injection Needle: SafetyGlide™ ¹ 25 Ga. x 1" (2.54 cm) | 1: Surgical Apparel: Bouffant Cap |
| 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm) | 1: Tubing±: 6" Extension with One-Way Stopcock |
| 1: Syringe: 3 mL Luer-Lock | ¹ A trademark of Becton, Dickinson and Company. |
| 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep | ² A registered trademark of CareFusion or one of its subsidiaries. |
| 1: 3 mL Applicator 2% CHG and 70% IPA Chloraprep® ² One-Step Solution with Hi-Lite Orange™ Tint | ³ A registered trademark of TIDI. |
| 1: SharpsAway® Disposal Cup | ⁴ A registered trademark of Johnson & Johnson Corporation. |
| 1: SharpsAway® II Locking Disposal Cup | |
| 1: Drape: 36" x 41" (91 cm x 104 cm) with 3" x 5" (8 cm x 13 cm) fenestration | |
| 1: Towel: 13" x 18" (33 cm x 46 cm) | |
| 1: Checklist Sheet | |
| 1: Medication Label: 1% Lidocaine | |
| 1: Dressing: GRIP-LOK® ³ Catheter Fixation Device | |
| 4: Gauze Pad: 2" x 2" (5 cm x 5 cm) | |
| 4: Gauze Pad: 4" x 4" (10 cm x 10 cm) | |
| 1: Dressing: I.V. Advanced 10 cm x 12 cm | |
| 1: Dressing: BIOPATCH® ⁴ | |

Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU

California Prop. 65

California Prop. 65
⚠ WARNING: Cancer and Reproductive Harm.
www.P65Warnings.ca.gov

±Precaution: Contains Phthalate: DEHP.

Warning: The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction. Proceeding from the present state of scientific knowledge, risks for male premature infants cannot be excluded in the case of long-term exposure or application.

Medical products containing phthalates should be used only temporarily with pregnant women, nursing mothers, babies and infants.

Not made with natural rubber latex.

Store between 20–25°C (68–77°F).

Store between 20 – 25°C (68 – 77°F).
Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.

Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.



LOT



 Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

LB 059317 R01 (2022-02)



20 ^{Ga.} **3.81** ^{cm}
catheter
length

.018 ^{inch dia.}
spring-wire
guide

ARROW