



REF AK-04510-S

20<sup>Ga.</sup>

12<sup>cm catheter length</sup>

.025<sup>inch dia. spring-wire guide</sup>


# Arterial Line

Contents:

- 1: Indwelling Catheter: 20 Ga. x 12 cm Radiopaque FEP
- 1: Spring-Wire Guide: .025" (0.64 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on One End - "J" Tip on Other)
- 1: Injection Needle: 25 Ga. x 1" (2.54 cm)
- 1: Injection Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Introducer Needle: 20 Ga. x 2" (5.08 cm) XTW
- 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 2: Syringe: 3 mL Luer-Lock
- 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep
- 1: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®<sup>1</sup> One-Step Solution with Hi-Lite Orange™ Tint
- 1: Drape: 18" x 26" (46 cm x 66 cm) with 3" (7.6 cm) fenestration
- 1: Scalpel: #11
- 3: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 1: Suture: 3-0 Braided Silk with Straight Needle

<sup>1</sup>A registered trademark of CareFusion or one of its subsidiaries.

**Rx only**  
**Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. [www.teleflex.com/IFU](http://www.teleflex.com/IFU)**  
**California Prop. 65**  
**⚠ WARNING:** Cancer and Reproductive Harm.  
[www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)  
**Not made with natural rubber latex.**  
Store between 20 – 25°C (68 – 77°F).  
Sterilized using ethylene oxide unless otherwise indicated in the contents list.  
Fluid path components are non-pyrogenic.



Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA  
**Teleflex®**  
[www.teleflex.com](http://www.teleflex.com)  
LBL059322 R02 (2022-06)



REF AK-04510-S

20<sup>Ga.</sup>

12<sup>cm catheter length</sup>

.025<sup>inch dia. spring-wire guide</sup>

# Arterial Line