



## **Femoral Arterial Line Catheterization**

## Contents:

- Indwelling Catheter: 18 Ga. x 16 cm with Integral Suture Wing Spring-Wire Guide: .025" (0.64 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End "." Tip on Other) Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) Introducer Needle: 26 Ga. x 1" (2.54 cm) Introducer Needle: 20 Ga. x 1-1/2" (3.81 cm) XTW Introducer Needle: 18 Ga. x 3" (7.62 cm) RW Blunt Fill Needle: 18 Ga. x 3" (7.62 cm) RW Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm) Syringe: 5 mL Luer-Slip 5 mL 1% Lidocaine HCI Solution and Alcohol Prep Drape: 24" x 36" (61 cm x 91 cm) with 4" (10 cm) fenestration Scalpel: #11 Prep Razor Prep Sponge Swab Gauze Pad: 4" x 4" (10 cm x 10 cm) Gauze Pad: 4" x 4" (10 cm x 10 cm) Suture: 3-0 Braided Silk with Straight Needle Tray: Prep 1: 1:

- 2:

- 1: 2: 4: 8: 1: 1:

Store between  $20 - 25^{\circ}C$  (68 - 77°F). Sterilized using ethylene oxide unless otherwise indicated in the contents list. Fluid path components are non-pyrogenic.

LOT

 $\sum$ 

Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA

**Neleflex**® LBL059324 R01 (2022-07)



~~~[

18<sup>Ga.</sup> 16<sup>cm</sup><sub>catheter</sub> .025<sup>inch dia.</sup> spring-wire

## **Femoral Arterial Line Catheterization**

