



REF ASK-00820-CHUS

20 Ga.

8 cm catheter length

.021 inch dia. spring-wire guide

Arterial Catheterization Kit

- Catheter: 20 Ga. x 8 cm Radiopaque Polyurethane
 Spring-Wire Guide, Marked: .021" (0.53 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on Both Ends)
 Introducer Needle: 20 Ga. x 1-1/2" (4 cm) XTW
 Injection Needle: SafetyGlide™1 18 Ga. x 1-1/2" (3.81 cm)
 Injection Needle: SafetyGlide™1 25 Ga. x 1" (2.54 cm)
 Syringe: 10 mL Luer-Lock
 Syringe: 5 mL Luer-Lock
 Syringe: 5 mL Luer-Slip
 Dust Cap: Non-Vented
 SharpsAway® Disposal Cup
 SharpsAway® Disposal Cup
 Drape: 41" x 55" with 4" fenestration
 Needle Holder
 Safety Scalpel: #11
 Gauze Pad: 4" x 4" (10 cm x 10 cm)
 Transducer Cover: 14 cm x 147 cm
 Suture: 3-0 Silk with Curved Needle
 Tubing±: 24" (61 cm) Pressure Monitoring

¹A trademark of Becton, Dickinson and Company

Catheter Priming Volume: 0.17 mL Priming volume is approximate.

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

California Prop. 65

Marning: Cancer and Reproductive Harm.

Www.P65Warnings.ca.gov

Precaution: Contains Phthalate: DEHP.

Warning: The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction. Proceeding from the present state of scientific knowledge, risks for male premature infants cannot be excluded in the case of long-term exposure or application. Medical products containing phthalates should be used only temporarily with pregnant women, nursing mothers, babies and infants.

Not made with natural rubber latex.

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.

Fluid path components are non-pyrogenic.



















Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA



LBL059335 R02 (2023-05)









