



ARROW®

REF AK-04020-C

20<sup>Ga.</sup> 4.45<sup>cm</sup> catheter length .018<sup>inch dia.</sup> spring-wire guides

## Radial Artery Catheterization

### Contents:

- 1: Catheter: 20 Ga. x 4.45 cm Radiopaque Polyurethane over 22 Ga. TW Introducer Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide
- 1: Injection Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Injection Needle: 27 Ga. x 1/2" (1.27 cm)
- 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Syringe: 3 mL Luer-Lock
- 1: Foil Pack 70% Alcohol Prep Pad
- 1: Skin Protectant Prep
- 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep
- 1: 1.5 mL Applicator 2% CHG and 70% IPA Chloraprep®<sup>1</sup> One-Step Frepp® Solution
- 1: Wing Clip: Suture
- 1: Towel: 13" x 18" (33 cm x 46 cm)
- 1: Drape: 18" x 26" with 2" x 3" fenestration, with adhesive
- 1: Scalpel: #11
- 4: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 4: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Suture: 3-0 Braided Silk with Straight Needle

<sup>1</sup>A registered trademark of CareFusion or one of its subsidiaries.

**Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

#### California Prop. 65

**⚠ WARNING:** Cancer and Reproductive Harm.  
[www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)

**Not made with natural rubber latex.**

Store between 20 – 25°C (68 – 77°F).

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.



ARROW®

REF AK-04020-C

20<sup>Ga.</sup> 4.45<sup>cm</sup> catheter length .018<sup>inch dia.</sup> spring-wire guides

## Radial Artery Catheterization

### Contents:

- 1: Catheter: 20 Ga. x 4.45 cm Radiopaque Polyurethane over 22 Ga. TW Introducer Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide
- 1: Injection Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Injection Needle: 27 Ga. x 1/2" (1.27 cm)
- 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Syringe: 3 mL Luer-Lock
- 1: Foil Pack 70% Alcohol Prep Pad
- 1: Skin Protectant Prep
- 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep
- 1: 1.5 mL Applicator 2% CHG and 70% IPA Chloraprep®<sup>1</sup> One-Step Frepp® Solution
- 1: Wing Clip: Suture
- 1: Towel: 13" x 18" (33 cm x 46 cm)
- 1: Drape: 18" x 26" with 2" x 3" fenestration, with adhesive
- 1: Scalpel: #11
- 4: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 4: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Suture: 3-0 Braided Silk with Straight Needle

<sup>1</sup>A registered trademark of CareFusion or one of its subsidiaries.

**Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

#### California Prop. 65

**⚠ WARNING:** Cancer and Reproductive Harm.  
[www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)

**Not made with natural rubber latex.**

Store between 20 – 25°C (68 – 77°F).

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.



LOT



Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA

**Teleflex®**  
[www.teleflex.com](http://www.teleflex.com)

LBL059341 R01 (2023-03)

LOT



Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA

**Teleflex®**  
[www.teleflex.com](http://www.teleflex.com)

LBL059341 R01 (2023-03)



REF AK-04020-C

ARROW®

20<sup>Ga.</sup> 4.45<sup>cm</sup> catheter length .018<sup>inch dia.</sup> spring-wire guides

## Radial Artery Catheterization



REF AK-04020-C

ARROW®

20<sup>Ga.</sup> 4.45<sup>cm</sup> catheter length .018<sup>inch dia.</sup> spring-wire guides

## Radial Artery Catheterization