



REF SAC-00520-X1A

20 Ga.

5 cm catheter length

.021 inch dia. spring-wire guide

Arterial Catheterization

Contents:

- 1: Catheter: 20 Ga. x 5 cm Radiopaque PEBA
- 1: Spring-Wire Guide, Marked: .021" (0.53 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on Both Ends)
- 1: Introducer Needle: 20 Ga. x 1-1/2" (4 cm) XTW
- 1: Injection Needle: SafetyGlide™ 1 25 Ga. x 1" (2.54 cm)
- 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Syringe: 3 mL Luer-Lock
- 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep
- 1: 3 mL Applicator 2% CHG and 70% IPA Chloraprep®² One-Step Solution with Hi-Lite Orange™ Tint
- 1: SharpsAway® Disposal Cup
- 1: SharpsAway® II Locking Disposal Cup
- 1: Drape: 36" x 41" (91 cm x 104 cm) with 3" x 5" (8 cm x 13 cm) fenestration
- 1: Towel: 13" x 18" (33 cm x 46 cm)
- 1: Needle Holder
- 1: Safety Scalpel: #11
- 1: Checklist Sheet
- 1: Medication Label: 1% Lidocaine
- 4: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 4: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Dressing: Tegaderm®³ 10 cm x 12 cm
- 1: Dressing: BIOPATCH®⁴
- 1: Surgical Apparel: Bouffant Cap
- 1: Surgical Apparel: Mask with Eye Shield
- 1: Suture: 3-0 Silk with Curved Needle

¹A trademark of Becton, Dickinson and Company.
²A registered trademark of CareFusion or one of its subsidiaries.
³A registered trademark of 3M Company.
⁴A registered trademark of Johnson & Johnson Corporation.

Catheter Priming Volume: 0.17 mL Priming volume is approximate.

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU

California Prop. 65
 ⚠ **WARNING:** Cancer and Reproductive Harm.
www.P65Warnings.ca.gov

Not made with natural rubber latex.
 Store between 20 – 25°C (68 – 77°F).
 Sterilized using ethylene oxide unless otherwise indicated in the contents list.
 Fluid path components are non-pyrogenic.
 Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.



Arrow International LLC
 Subsidiary of Teleflex Incorporated
 3015 Carrington Mill Blvd.
 Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

LBL059367 R00 (2022-01)



REF SAC-00520-X1A

20 Ga.

5 cm catheter length

.021 inch dia. spring-wire guide

Arterial Catheterization