



REF NA-04020-X1ASP

20 Ga. 4.45 cm catheter length .018 inch dia. spring-wire guide

Radial Artery Catheterization with Integral Needle Protection

Contents:

- | | |
|--|---|
| 1: Catheter: 20 Ga. x 4.45 cm Radiopaque Polyurethane over 22 Ga. TW Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide | 1: Surgical Apparel: Mask with Eye Shield |
| 1: Injection Needle: SafetyGlide™ 1 25 Ga. x 1" (2.54 cm) | 1: Suture: 3-0 Silk with Curved Needle |
| 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm) | 1: Tubing±: 6" Extension with One-Way Stopcock |
| 1: Syringe: 3 mL Luer-Lock | ¹ A trademark of Becton, Dickinson and Company. |
| 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep | ² A registered trademark of CareFusion or one of its subsidiaries. |
| 1: 3 mL Applicator 2% CHG and 70% IPA Chloraprep® ² One-Step Solution with Hi-Lite Orange™ Tint | ³ A registered trademark of 3M Company. |
| 1: SharpsAway® Disposal Cup | ⁴ A registered trademark of Johnson & Johnson Corporation. |
| 1: SharpsAway® II Locking Disposal Cup | |
| 1: Drape: 36" x 41" (91 cm x 104 cm) with 3" x 5" (8 cm x 13 cm) fenestration | |
| 1: Towel: 13" x 18" (33 cm x 46 cm) | |
| 1: Needle Holder | |
| 1: Safety Scalpel: #11 | |
| 1: Checklist Sheet | |
| 1: Medication Label: 1% Lidocaine | |
| 4: Gauze Pad: 2" x 2" (5 cm x 5 cm) | |
| 4: Gauze Pad: 4" x 4" (10 cm x 10 cm) | |
| 1: Dressing: Tegaderm® ³ 10 cm x 12 cm | |
| 1: Dressing: BIOPATCH® ⁴ | |
| 1: Surgical Apparel: Bouffant Cap | |

Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU

California Prop. 65

⚠ **WARNING:** Cancer and Reproductive Harm. www.P65Warnings.ca.gov

±**Precaution: Contains Phthalate: DEHP.**

Warning: The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction. Proceeding from the present state of scientific knowledge, risks for male premature infants cannot be excluded in the case of long-term exposure or application. Medical products containing phthalates should be used only temporarily with pregnant women, nursing mothers, babies and infants.


Not made with natural rubber latex.

Store between 20 – 25°C (68 – 77°F).

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.

Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.



Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

LBL059374 R00 (2021-11)



REF NA-04020-X1ASP

20 Ga. 4.45 cm catheter length .018 inch dia. spring-wire guide

Radial Artery Catheterization with Integral Needle Protection