



REF SAC-00522-1A

22 Ga.

5 cm catheter length

.021 inch dia. spring-wire guide

# Arterial Catheterization

**Contents:**

- |  |   |
|--|---|
| <ul style="list-style-type: none"> <li>1: Catheter: 22 Ga. x 5 cm Radiopaque PEBA</li> <li>1: Spring-Wire Guide, Marked: .021" (0.53 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on Both Ends)</li> <li>1: Introducer Needle: 22 Ga. x 1-1/2" (4 cm) XTW</li> <li>1: Injection Needle: SafetyGlide™ 25 Ga. x 1" (2.54 cm)</li> <li>1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)</li> <li>1: Syringe: 3 mL Luer-Lock</li> <li>1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep</li> <li>1: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®<sup>2</sup> One-Step Solution with Hi-Lite Orange™ Tint</li> <li>1: SharpsAway® Disposal Cup</li> <li>1: SharpsAway® II Locking Disposal Cup</li> <li>1: Towel: 13" x 18" (33 cm x 46 cm)</li> <li>1: Drape: 18" x 26" with 2" x 3" fenestration, with adhesive</li> <li>1: Needle Holder</li> <li>1: Safety Scalpel: #11</li> <li>1: Medication Label: 1% Lidocaine</li> <li>4: Gauze Pad: 2" x 2" (5 cm x 5 cm)</li> <li>4: Gauze Pad: 4" x 4" (10 cm x 10 cm)</li> <li>1: Dressing: Tegaderm®<sup>3</sup> 10 cm x 12 cm</li> <li>1: Surgical Apparel: Bouffant Cap</li> </ul> | <ul style="list-style-type: none"> <li>1: Surgical Apparel: Mask with Eye Shield</li> <li>1: Suture: 3-0 Silk with Curved Needle</li> </ul> <p><sup>1</sup>A trademark of Becton, Dickinson and Company.<br/> <sup>2</sup>A registered trademark of CareFusion or one of its subsidiaries.<br/> <sup>3</sup>A registered trademark of 3M Company.</p> |
|--|---|

**Catheter Priming Volume: 0.16 mL Priming volume is approximate.**

**Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. [www.teleflex.com/IFU](http://www.teleflex.com/IFU)**  
**California Prop. 65**

**⚠ WARNING:** Cancer and Reproductive Harm.  
[www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)

**Not made with natural rubber latex.**

Store between 20 – 25°C (68 – 77°F).

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.

Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.



Arrow International LLC  
 Subsidiary of Teleflex Incorporated  
 3015 Carrington Mill Blvd.  
 Morrisville, NC 27560 USA

**Teleflex®**  
[www.teleflex.com](http://www.teleflex.com)

LBL059378 R00 (2022-01)



REF SAC-00522-1A

22 Ga.

5 cm catheter length

.021 inch dia. spring-wire guide

# Arterial Catheterization