



REF SAC-00522-X1A

22^{Ga.}

5^{cm}
catheter
length

.021^{inch dia.}
spring-wire
guide

Arterial Catheterization

Contents:

- | | |
|---|--|
| <ul style="list-style-type: none"> 1: Catheter: 22 Ga. x 5 cm Radiopaque PEBA 1: Spring-Wire Guide, Marked: .021" (0.53 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on Both Ends) 1: Introducer Needle: 22 Ga. x 1-1/2" (4 cm) XTW 1: Injection Needle: SafetyGlide™¹ 25 Ga. x 1" (2.54 cm) 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm) 1: Syringe: 3 mL Luer-Lock 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep 1: 3 mL Applicator 2% CHG and 70% IPA Chloraprep®² One-Step Solution with Hi-Lite Orange™ Tint 1: SharpsAway® Disposal Cup 1: SharpsAway® II Locking Disposal Cup 1: Towel: 13" x 18" (33 cm x 46 cm) 1: Drape: 18" x 26" with 2" x 3" fenestration, with adhesive 1: Needle Holder 1: Safety Scalpel: #11 1: Checklist Sheet 1: Medication Label: 1% Lidocaine 4: Gauze Pad: 2" x 2" (5 cm x 5 cm) 4: Gauze Pad: 4" x 4" (10 cm x 10 cm) 1: Dressing: Tegaderm®³ 10 cm x 12 cm | <ul style="list-style-type: none"> 1: Dressing: BIOPATCH®⁴ 1: Surgical Apparel: Bouffant Cap 1: Surgical Apparel: Mask with Eye Shield 1: Suture: 3-0 Silk with Curved Needle <p>¹A trademark of Becton, Dickinson and Company.
 ²A registered trademark of CareFusion or one of its subsidiaries.
 ³A registered trademark of 3M Company.
 ⁴A registered trademark of Johnson & Johnson Corporation.</p> |
|---|--|

Catheter Priming Volume: 0.16 mL Priming volume is approximate.

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU
California Prop. 65

⚠ **WARNING:** Cancer and Reproductive Harm.
www.P65Warnings.ca.gov

Not made with natural rubber latex.

Store between 20 – 25°C (68 – 77°F).

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.

Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.



LOT



REF SAC-00522-X1A

22^{Ga.}

5^{cm}
catheter
length

.021^{inch dia.}
spring-wire
guide

Arterial Catheterization

Arrow International LLC
 Subsidiary of Teleflex Incorporated
 3015 Carrington Mill Blvd.
 Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

LBL059379 R00 (2022-01)