



REF SAC-00820-X1A (IPNXXXXXX)

20 Ga. 8 cm catheter length .021 inch dia. spring-wire guide

# Arterial Catheterization

### Contents:

- 1: Catheter: 20 Ga. x 8 cm Radiopaque PEBA Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX
- 1: Spring-Wire Guide, Marked: .021" (0.53 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on Both Ends) Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX
- 1: Introducer Needle: 20 Ga. x 1-1/2" (4 cm) XTW Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX
- 1: Injection Needle: SafetyGlide™ 25 Ga. x 1" (2.54 cm) Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX
- 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm) Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX
- 1: Syringe: 3 mL Luer-Lock Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX
- 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX
- 1: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®<sup>2</sup> One-Step Solution with Hi-Lite Orange™ Tint Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX
- 1: SharpsAway® Disposal Cup Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX
- 1: SharpsAway® II Locking Disposal Cup Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX
- 1: Drape: 36" x 41" (91 cm x 104 cm) with 3" x 5" (8 cm x 13 cm) fenestration Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX
- 1: Towel: 13" x 18" (33 cm x 46 cm) Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX
- 1: Needle Holder Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX
- 1: Safety Scalpel: #11 Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX
- 1: Checklist Sheet Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX
- 1: Medication Label: 1% Lidocaine Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX
- 4: Gauze Pad: 2" x 2" (5 cm x 5 cm) Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX
- 4: Gauze Pad: 4" x 4" (10 cm x 10 cm) Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX
- 1: Dressing: Tegaderm®<sup>3</sup> 10 cm x 12 cm Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX
- 1: Dressing: BIOPATCH®<sup>4</sup> Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX
- 1: Surgical Apparel: Bouffant Cap Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX
- 1: Surgical Apparel: Mask with Eye Shield Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX
- 1: Suture: 3-0 Silk with Curved Needle Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX

<sup>1</sup>A trademark of Becton, Dickinson and Company.  
<sup>2</sup>A registered trademark of CareFusion or one of its subsidiaries.  
<sup>3</sup>A registered trademark of 3M Company.  
<sup>4</sup>A registered trademark of Johnson & Johnson Corporation.

Catheter Priming Volume: 0.17 mL Priming volume is approximate.

**Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU California Prop. 65**

**⚠ WARNING:** Cancer and Reproductive Harm. www.P65Warnings.ca.gov

**Not made with natural rubber latex.**

Store between 20 – 25°C (68 – 77°F). Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic. Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.

Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA

Set Product of xxx  
Packaged in xxx

**Teleflex®**  
www.teleflex.com  
LBL059380 R01 (2025-05)



LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXXXXX  
(17)YYMMDD  
(11)YYMMDD  
(10)Lot Number



REF SAC-00820-X1A

20 Ga. 8 cm catheter length .021 inch dia. spring-wire guide



# Arterial Catheterization