



REF NA-04550-X1A

(IPNXXXXXX)

20<sup>Ga.</sup>12<sup>cm</sup>  
catheter  
length.025<sup>inch dia.</sup>  
spring-wire  
guide

ARROW  
ErgoPack®



# Arterial Catheterization

## Contents:

- |  |   |
|--|---|
| 1: One-Lumen Catheter: 20 Ga. (1.2 mm OD) x 12 cm  | 4: Gauze Pad: 4" x 4" (10 cm x 10 cm)   |
| 1: Spring-Wire Guide: .025" (0.64 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer | 1: Dressing: Tegaderm® <sup>3</sup> 10 cm x 12 cm                             |
| 1: Introducer Needle: 20 Ga. x 2" (5.08 cm) XTW  | 1: Surgical Apparel: Impervious Gown  |
| 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) RW   | 1: Dressing: BIOPATCH® <sup>4</sup>   |
| 1: Injection Needle: SafetyGlide™ <sup>1</sup> 25 Ga. x 1" (2.54 cm)   | 1: Surgical Apparel: Bouffant Cap   |
| 1: Injection Needle: SafetyGlide™ <sup>1</sup> 22 Ga. x 1-1/2" (3.81 cm) RW  | 1: Surgical Apparel: Mask with Eye Shield                                     |
| 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)  | 1: Suture: 3-0 Silk with Curved Needle  |
| 1: Syringe: 3 mL Luer-Lock   | 1: Tubing: 6" Extension with One-Way Stopcock                                 |
| 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep   | 1: HemoHopper® Fluid Receptacle   |
| 2: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep® <sup>2</sup> One-Step Solution with Hi-Lite Orange™ Tint                         | <sup>1</sup> A trademark of Becton, Dickinson and Company.                    |
| 1: SharpsAway® Disposal Cup  | <sup>2</sup> A registered trademark of CareFusion or one of its subsidiaries. |
| 1: SharpsAway® II Locking Disposal Cup   | <sup>3</sup> A registered trademark of 3M Company.                            |
| 1: Maximal Barrier Drape™ with 4" fenestration   | <sup>4</sup> A registered trademark of Johnson & Johnson Corporation.         |
| 1: Towel: 13" x 18" (33 cm x 46 cm)  |   |
| 1: Needle Holder   |   |
| 1: Safety Scalpel: #11   |   |
| 1: Checklist Sheet   |   |
| 1: Medication Label: 1% Lidocaine  |   |
| 4: Gauze Pad: 2" x 2" (5 cm x 5 cm)  |   |

±Precaution: Contains Phthalate: DEHP.

**Warning:** The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction. Proceeding from the present state of scientific knowledge, risks for male premature infants cannot be excluded in the case of long-term exposure or application. Medical products containing phthalates should be used only temporarily with pregnant women, nursing mothers, babies and infants.

**Not made with natural rubber latex.**

Store between 20 – 25°C (68 – 77°F).

Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.

