



REF EDC-00820-B

20 Ga. catheter

8 cm catheter length

.010 inch dia. wire guide

ARROW ENDURANCE™ Extended Dwell Peripheral Catheter System

Contents:

- 1: Catheter System: 20 Ga. (1.10 mm) x 80 mm (3.15") ARROW ENDURANCE™ over 22 Ga. (0.8 mm) RW Introducer Needle with Integral 0.010" (0.26 mm) Straight Soft Tip Wire Guide
- 1: Skin Protectant Prep
- 1: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®¹ One-Step Solution with Hi-Lite Orange™ Tint
- 1: Towel: 24" x 36" (61 cm x 91 cm)
- 1: Arm Drape: 36" x 46" with 3" x 5" Fenestration
- 1: Label: Venous/Arterial Placement
- 1: Paper Tape Measure
- 1: Tourniquet
- 4: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 4: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Dressing: Tegaderm®² 10 cm x 12 cm
- 1: Dressing: STATLOCK®³ Catheter Stabilization Device
- 1: Surgical Apparel: Bouffant Cap
- 1: Surgical Apparel: Mask

¹A registered trademark of CareFusion or one of its subsidiaries.

²A registered trademark of 3M Company.

³A registered trademark of C. R. Bard, Inc.

Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU

California Prop. 65


⚠ WARNING: Cancer and Reproductive Harm.

www.P65Warnings.ca.gov

Not made with natural rubber latex.

Store between 15 – 30°C (59 – 86°F). Avoid freezing and excessive heat above 40°C (104°F).

Fluid path components are non-pyrogenic.



LOT



Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

LBL060977 R00 (2022-02)



REF EDC-00820-B

20 Ga. catheter

8 cm catheter length

.010 inch dia. wire guide

ARROW ENDURANCE™ Extended Dwell Peripheral Catheter System