



REF EDC-00822-B

22 Ga. catheter

8 cm catheter length

.010 inch dia. wire guide

# ARROW ENDURANCE™ Extended Dwell Peripheral Catheter System

**Contents:**

- 1: Catheter System: 22 Ga. (0.91 mm) x 80 mm (3.15") ARROW ENDURANCE™ over 24 Ga. (0.60 mm) TW Introducer Needle with Integral 0.010" (0.26 mm) Straight Soft Tip Wire Guide
- 1: Skin Protectant Prep
- 1: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®<sup>1</sup> One-Step Solution with Hi-Lite Orange™ Tint
- 1: Towel: 24" x 36" (61 cm x 91 cm)
- 1: Arm Drape: 36" x 46" with 3" x 5" Fenestration
- 1: Label: Venous/Arterial Placement
- 1: Paper Tape Measure
- 1: Tourniquet
- 4: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 4: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Dressing: Tegaderm®<sup>2</sup> 10 cm x 12 cm
- 1: Dressing: STATLOCK®<sup>3</sup> Catheter Stabilization Device
- 1: Surgical Apparel: Bouffant Cap
- 1: Surgical Apparel: Mask

<sup>1</sup>A registered trademark of CareFusion or one of its subsidiaries.

<sup>2</sup>A registered trademark of 3M Company.

<sup>3</sup>A registered trademark of C. R. Bard, Inc.

**Rx only**

**Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. [www.teleflex.com/IFU](http://www.teleflex.com/IFU)**

**California Prop. 65**

**⚠ WARNING:** Cancer and Reproductive Harm.

[www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)

**Not made with natural rubber latex.**

Store between 15 – 30°C (59 – 86°F). Avoid freezing and excessive heat above 40°C (104°F).

Fluid path components are non-pyrogenic.



LOT



Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA

**Teleflex®**  
[www.teleflex.com](http://www.teleflex.com)

LBL060978 R00 (2022-02)



REF EDC-00822-B

22 Ga. catheter

8 cm catheter length

.010 inch dia. wire guide

# ARROW ENDURANCE™ Extended Dwell Peripheral Catheter System