



REF AK-09000 (IPNXXXXXX)

Peritoneal Lavage

8 Fr. 24 cm catheter length .035 inch dia. spring-wire guide

Contents:

- 1: Peritoneal Lavage Catheter: 8 Fr. x 9-1/4" (24 cm) Radiopaque FEP with Side Holes origin:xxxxxxxxxxxxxxxxxxxxxxxx
- 1: Spring-Wire Guide: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) origin:xxxxxxxxxxxxxxxxxxxxxxxx
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle origin:xxxxxxxxxxxxxxxxxxxxxxxx
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW origin:xxxxxxxxxxxxxxxxxxxxxxxx
- 1: Injection Needle: SafetyGlide™1 25 Ga. x 1" (2.54 cm) origin:xxxxxxxxxxxxxxxxxxxxxxxx
- 1: Injection Needle: SafetyGlide™1 22 Ga. x 1-1/2" (3.81 cm) RW origin:xxxxxxxxxxxxxxxxxxxxxxxx
- 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm) origin:xxxxxxxxxxxxxxxxxxxxxxxx
- 2: Syringe: 5 mL Luer-Slip origin:xxxxxxxxxxxxxxxxxxxxxxxx
- 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep origin:xxxxxxxxxxxxxxxxxxxxxxxx
- 1: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®2 One-Step Solution with Hi-Lite Orange™ Tint origin:xxxxxxxxxxxxxxxxxxxxxxxx
- 1: SharpsAway® II Locking Disposal Cup origin:xxxxxxxxxxxxxxxxxxxxxxxx
- 1: Drape: 18" x 26" (46 cm x 66 cm) with 3" (7.6 cm) fenestration origin:xxxxxxxxxxxxxxxxxxxxxxxx
- 1: Safety Scalpel: #11 origin:xxxxxxxxxxxxxxxxxxxxxxxx
- 1: Medication Label: 1% Lidocaine origin:xxxxxxxxxxxxxxxxxxxxxxxx
- 2: Gauze Pad: 2" x 2" (5 cm x 5 cm) origin:xxxxxxxxxxxxxxxxxxxxxxxx
- 4: Gauze Pad: 4" x 4" (10 cm x 10 cm) origin:xxxxxxxxxxxxxxxxxxxxxxxx
- 1: Tubing Set: Fluid Administration origin:xxxxxxxxxxxxxxxxxxxxxxxx

¹A trademark of Becton, Dickinson and Company.
²A registered trademark of CareFusion or one of its subsidiaries.

Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU

California Prop. 65

⚠ WARNING: Cancer and Reproductive Harm.
www.P65Warnings.ca.gov

±Precaution: Contains Phthalate: DEHP.

Warning: The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction. Proceeding from the present state of scientific knowledge, risks for male premature infants cannot be excluded in the case of long-term exposure or application. Medical products containing phthalates should be used only temporarily with pregnant women, nursing mothers, babies and infants.

Not made with natural rubber latex.

Store between 20 – 25°C (68 – 77°F).

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.

Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.

Contraindications: Patients with an acute abdomen requiring immediate surgery or penetrating abdominal trauma.



LOT LotNumber

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)LotNumber

Set Product of xxx
Packaged in xxx

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA



LBL060981 R01 (2024-04)

REF AK-09000

8 Fr. 24 cm catheter length .035 inch dia. spring-wire guide

Peritoneal Lavage

