

REF AK-09001

# Peritoneal Lavage

8 Fr.

24 cm catheter length

.035 inch dia. spring-wire guide

**Contents:**

- 1: Peritoneal Lavage Catheter: 8 Fr. x 9-1/4" (24 cm) Radiopaque FEP with Side Holes
- 1: Spring-Wire Guide: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other)
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)
- 1: Injection Needle: 25 Ga. x 1" (2.54 cm)
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 2: Syringe: 5 mL Luer-Slip
- 1: Drape: 18" x 26" (46 cm x 66 cm) with 3" (7.6 cm) fenestration
- 1: Scalpel: #11
- 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 4: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Tubing Set: Fluid Administration

**Rx only**

**Warning:** Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. [www.teleflex.com/IFU](http://www.teleflex.com/IFU)

California Prop. 65

**WARNING:** Cancer and Reproductive Harm. [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)

**Precaution:** Contains Phthalate: DEHP.

**Warning:** The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction. Proceeding from the present state of scientific knowledge, risks for male premature infants cannot be excluded in the case of long-term exposure or application. Medical products containing phthalates should be used only temporarily with pregnant women, nursing mothers, babies and infants.

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.

**Contraindications:** Patients with an acute abdomen requiring immediate surgery or penetrating abdominal trauma.



LOT



Arrow International LLC  
 Subsidiary of Teleflex Incorporated  
 3015 Carrington Mill Blvd.  
 Morrisville, NC 27560 USA

**Teleflex**  
[www.teleflex.com](http://www.teleflex.com)

LBL061013 R00 (2022-03)

REF AK-09001

8 Fr.

24 cm catheter length

.035 inch dia. spring-wire guide

# Peritoneal Lavage