



REF ASK-09903-SM1 (IPNXXXXXX)

ASK ARROW®

9^{Fr.}10^{cm}
sheath
length.035^{inch dia.}
spring-wire
guide

Percutaneous Sheath Introducer Procedural Kits With Silk Sutures



Contents:

- | | |
|---|--|
| 1: Arrow-Flex® Sheath: 9 Fr. (2.9 mm tip ID) x 10 cm and Tissue Dilator | 1: Towel |
| 1: Spring-Wire Guide, Nitinol, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow® GlideWheel™ Wire Advancer | 1: Needle Holder |
| 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle | 1: Forceps |
| 1: Pressure Transduction Probe | 1: Safety Scalpel: #11 |
| 1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe | 1: Stopcock: 4-way High-Flow |
| 1: Injection Needle: SafetyGlide™ ¹ 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock Syringe | 1: Checklist/CLIP Sheet |
| 1: Injection Needle: SafetyGlide™ ¹ 23 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe | 1: Sterile Procedure Sign |
| 1: Syringe: 10 mL Luer-Lock | 2: Gauze Pad: 2" x 2" (5 cm x 5 cm) |
| 1: Obturator: 8 Fr. | 5: Gauze Pad: 4" x 4" (10 cm x 10 cm) |
| 1: Cath-Gard® Catheter Contamination Shield: 80 cm with TwistLock™ Adapter | 1: Dressing: I.V. Advanced 10 cm x 12 cm |
| 1: Dust Cap: Non-Vented | 1: Transducer Cover: 14 cm x 147 cm |
| 1: SharpsAway® II Locking Disposal Cup | 1: Suture: 3-0 Braided Silk with Straight Needle |
| 1: SharpsAway® Disposal Cup | 1: Suture: 3-0 Silk with Curved Needle |
| 1: Maximal Barrier Drape™ with 4" fenestration | 1: Tubing±: 20" (50 cm) Extension |
| 1: Drape: 21" x 36" (53 cm x 91 cm) with 3" (7.6 cm) fenestration, with adhesive | 1: HemoHopper® Fluid Receptacle |

¹A trademark of Becton, Dickinson and Company.

±Precaution: Contains Phthalate: DEHP.

Warning: The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction. Proceeding from the present state of scientific knowledge, risks for male premature infants cannot be excluded in the case of long-term exposure or application. Medical products containing phthalates should be used only temporarily with pregnant women, nursing mothers, babies and infants.

Not made with natural rubber latex.

