


# ARROW™ EZ-IO™

INTRAOSSSEOUS VASCULAR ACCESS



## Intraosseous Access Tray + Single-Use Power Driver

**REF 9018T-VC-005**  
(IPNXXXXXX)

**15 mm 15 ga** **15 mm x 1.8 mm**  
3-39 kg

**Contents:**

- 1: EZ-IO™ Power Driver (Single-Use)
  - 1: EZ-IO™ 15mm 15 Ga Intraosseous Needle‡ Set
  - 1: 10 mL Luer-Lock Pre-Filled Saline Syringe
  - 1: 1.5 mL Applicator 2% CHG and 70% IPA ChloroPrep®<sup>1</sup> One-Step Frepp® Solution
  - 1: SharpsAway® II Locking Disposal Cup
  - 1: EZ-IO™ Wristband
  - 1: EZ-Stabilizer™ Dressing
  - 1: EZ-Connect™ Extension Set (Priming Volume - Approx. 1.0 ml)
- <sup>1</sup>A registered trademark of CareFusion or one of its subsidiaries.

Teleflex Medical  
IDA Business and  
Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

**Teleflex™**  
ArrowEZIO.com  
24-Hour Clinical Support: 1.888.413.3104  
Customer Service: 1.866.479.8500

Product of xxx  
Packaged in xxx

**Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

**California Prop. 65**  
⚠ **WARNING:** Cancer and Reproductive Harm.  
[www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)

**‡MR Conditional**  
**Not made with natural rubber latex.**  
Store between 15 – 30°C (59 – 86°F). Avoid freezing and excessive heat above 40°C (104°F).

Rx only

**LOT** Lot Number  
 YYYY-MM-DD YYYY-MM-DD

(01)XXXXXXXXXXXXXX  
(17)YYMMDD  
(11)YYMMDD  
(10)Lot Number

LBL069290 R02 (2024-01)

Peel Here

**ARROW™ EZ-IO™**  
INTRAOSSSEOUS VASCULAR ACCESS

**REF 9018T-VC-005**

**Intraosseous Access Tray**  
+ Single-Use Power Driver

**15 mm 15 ga** **15 mm x 1.8 mm**  
3-39 kg