



REF **CDA-26142-XCN1A** (IPNXXXXXX)

**14** Fr. **2** Lumen **25** cm catheter length **.035** inch dia. spring-wire guide



# Arrowg+ard Blue® Two-Lumen Hemodialysis Catheterization for High Volume Infusions



**Contents:**

- |   |   |
|---|---|
| <ul style="list-style-type: none"> <li>1: Two-Lumen Catheter: 14 Fr. (4.72 mm OD) x 25 cm Radiopaque Polyurethane with Blue FlexTip®, Arrowg+ard® Antimicrobial Surface Treatment, Extension Line Clamps</li> <li>1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow® GlideWheel™ Wire Advancer</li> <li>1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle</li> <li>1: Pressure Transduction Probe</li> <li>1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe</li> <li>1: Injection Needle: Eclipse™<sup>1</sup> 23 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe</li> <li>1: Injection Needle: Eclipse™<sup>1</sup> 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock Syringe</li> <li>1: Syringe: 10 mL Luer-Lock</li> <li>1: Tissue Dilator: 14 Fr. (4.7 mm) x 14 cm</li> <li>1: Tissue Dilator: 12 Fr. (4.0 mm) x 14 cm</li> <li>1: SharpsAway® II Locking Disposal Cup</li> <li>1: SharpsAway® Disposal Cup</li> <li>1: Maximal Barrier Drape™ with 4" fenestration</li> <li>1: Towel</li> </ul> | <ul style="list-style-type: none"> <li>1: Needle Holder</li> <li>1: Forceps</li> <li>1: Safety Scalpel: #11</li> <li>1: Checklist/CLIP Sheet</li> <li>1: Sterile Procedure Sign</li> <li>2: Gauze Pad: 2" x 2" (5 cm x 5 cm)</li> <li>5: Gauze Pad: 4" x 4" (10 cm x 10 cm)</li> <li>1: Surgical Apparel: Impervious Gown</li> <li>1: Dressing: Tegaderm®<sup>2</sup> CHG 3-1/2" x 4-1/2" (8.5 cm x 11.5 cm)</li> <li>1: Transducer Cover: 14 cm x 147 cm</li> <li>1: Surgical Apparel: Mask with Eye Shield</li> <li>1: Surgical Apparel: Bouffant Cap</li> <li>1: Suture: 3-0 Silk with Curved Needle</li> <li>1: HemoHopper® Fluid Receptacle</li> </ul> <p><sup>1</sup>A trademark of Becton, Dickinson and Company.<br/> <sup>2</sup>A registered trademark of 3M Company.</p> |
|---|---|

**Not made with natural rubber latex.**

Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F).

**Contraindications:** The Arrow large-bore two-lumen catheter is not designed for long-term (≥ 30 days) hemodialysis or for use in patients with thrombosed vessels.


The Arrowg+ard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine, and/or sulfa drugs.



**Priming Volume\***

Distal	1.80 mL
Proximal	1.60 mL

\*Priming volumes are approximate and are done without accessories.



† In vitro tests performed using 3:1 glycerine/saline at 37°C.

