



REF CDC-29903-1A

9 Fr.

10 cm sheath length

.035 inch dia. spring-wire guide

ARROW

ErgoPack

# Arrowgard Blue<sup>®</sup> PSI for use with 7.5 - 8 Fr. Catheters

## Contents:

- |  |  |
|--|--|
| 1: Arrowgard Blue <sup>®</sup> Arrow-Flex <sup>®</sup> Sheath: 9 Fr. (2.9 mm tip ID) x 10 cm and Tissue Dilator                            | 1: Maximal Barrier Drape <sup>™</sup> with 4" fenestration |
| 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer | 1: Needle Holder   |
| 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle   | 1: Safety Scalpel: #11                                     |
| 1: Injection Needle: SafetyGlide <sup>™</sup> 1 25 Ga. x 1" (2.54 cm)  | 1: Stopcock: 4-way High-Flow                               |
| 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)  | 1: Checklist/CLIP Sheet                                    |
| 1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Luer-Slip Syringe   | 1: Sterile Procedure Sign                                  |
| 1: Injection Needle: SafetyGlide <sup>™</sup> 1 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe                                       | 1: Medication Label: 1% Lidocaine                          |
| 1: Syringe: 3 mL Luer-Lock   | 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)                        |
| 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep   | 5: Gauze Pad: 4" x 4" (10 cm x 10 cm)                      |
| 1: 3 mL Applicator 2% CHG and 70% IPA Chloraprep <sup>®2</sup> One-Step Solution with Hi-Lite Orange <sup>™</sup> Tint                     | 1: Surgical Apparel: Impervious Gown                       |
| 1: Obturator: 8 Fr.  | 1: Surgical Apparel: Bouffant Cap                          |
| 1: Cath-Gard <sup>®</sup> Catheter Contamination Shield: 80 cm with TwistLock <sup>™</sup> Adapter   | 1: Surgical Apparel: Mask with Eye Shield                  |
| 1: Dust Cap: Non-Vented  | 1: Suture: 3-0 Silk with Curved Needle                     |
| 1: SharpsAway <sup>®</sup> II Locking Disposal Cup   | 1: Tubing: 20" (50 cm) Extension                           |
| 1: SharpsAway <sup>®</sup> Disposal Cup  | 1: HemoHopper <sup>®</sup> Fluid Receptacle                |
- <sup>1</sup>A trademark of Becton, Dickinson and Company.  
<sup>2</sup>A registered trademark of CareFusion or one of its subsidiaries.

±Precaution: Contains Phthalate: DEHP.

**Warning:** The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction. Proceeding from the present state of scientific knowledge, risks for male premature infants cannot be excluded in the case of long-term exposure or application. Medical products containing phthalates should be used only temporarily with pregnant women, nursing mothers, babies and infants.

Not made with natural rubber latex.

Store between 20 – 25°C (68 – 77°F).

Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.

**Contraindications:** The Arrowgard Blue antimicrobial sheath introducer is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine, and/or sulfa drugs (refer to product instructions for references).

**This Product Does Not Contain a 5 mL Arrow<sup>®</sup> Raulerson Syringe or Pressure Transduction Probe Due to a Supply Shortage, and Contains An Additional 5 mL Syringe.**

