



REF ASK-29903-HGH1

(IPNXXXXXX)

9 Fr.

10 cm sheath length

.035 inch dia. spring-wire guide



# Arrowgard Blue® PSI Kit for use with 7.5 - 8 Fr. Catheters

**Contents:**

- |  |   |
|--|---|
| 1: Arrowgard Blue® Arrow-Flex® Sheath: 9 Fr. (2.9 mm tip ID) x 10 cm and Tissue Dilator  | 5: Gauze Pad: 4" x 4" (10 cm x 10 cm)                               |
| 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer | 1: Dressing: Tegaderm® <sup>2</sup> CHG 4" x 4-3/4" (10 cm x 12 cm) |
| 1: Injection Needle: SafetyGlide™ <sup>1</sup> 25 Ga. x 1" (2.54 cm)   | 1: Suture: 3-0 Silk with Curved Needle                              |
| 1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe                   | 1: HemoHopper® Fluid Receptacle                                     |
| 1: Injection Needle: SafetyGlide™ <sup>1</sup> 23 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe  | <sup>1</sup> A trademark of Becton, Dickinson and Company.          |
| 1: Obturator: 8 Fr.  | <sup>2</sup> A registered trademark of 3M Company.                  |
| 1: Cath-Gard® Catheter Contamination Shield: 80 cm with TwistLock™ Adapter   |   |
| 1: Dust Cap: Non-Vented  |   |
| 1: SharpsAway® II Locking Disposal Cup   |   |
| 1: SharpsAway® Disposal Cup  |   |
| 1: Maximal Barrier Drape™ with 4" fenestration   |   |
| 1: Needle Holder   |   |
| 1: Safety Scalpel: #11   |   |
| 1: Stopcock: 4-way High-Flow   |   |
| 1: Checklist/CLIP Sheet  |   |
| 1: Sterile Procedure Sign  |   |
| 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)  |   |

**Not made with natural rubber latex.**

Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F).

**Contraindications:** The Arrowgard Blue antimicrobial sheath introducer is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine, and/or sulfa drugs (refer to product instructions for references).