

DO NOT BEND

ARROW®



REF SAC-00818-PBX

(IPNXXXXXX)

18^{Ga.}8^{cm}
catheter
length

Arterial Catheterization Set

Contents:

- 1: Catheter: 18 Ga. x 8 cm Radiopaque PEBA Origin: XX XXXX XX XXX XXXXXXXX
XXXXXXXXXX
- 1: Spring-Wire Guide, Marked: .025" (0.64 mm) dia. x 13-1/8" (33 cm) (Straight Soft Tip
on Both Ends) Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX
- 1: Introducer Needle: 18 Ga. x 2" (5 cm) XTW Origin: XX XXXX XX XXX XXXXXXXX
XXXXXXXXXX

All components are CE 2797 unless otherwise noted.

Catheter Priming Volume: 0.19 mL Priming volume is approximate.

Rx only Warning: Read all package insert warnings, precautions, and instructions
prior to use. Failure to do so may result in severe patient injury or death.

www.teleflex.com/IFU

California Prop. 65

⚠ **WARNING:** Cancer and Reproductive Harm.

www.P65Warnings.ca.gov

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.



CAS # 7440-48-4

EU Authorized
Representative and
Importer:

EC REP



Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

All components are CE 2797 unless otherwise noted.

LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD

Set Product of xxx
Packaged in xxx

(01)XXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number

Teleflex®
www.teleflex.com

LBL075263 R00 (2025-05)



REF SAC-00818-PBX

18^{Ga.}8^{cm}
catheter
length

Arterial Catheterization Set

ARROW®

DO NOT BEND