(IPNXXXXXX)





REF ASK-09907-SB1

10 sheath length

.025 inch dia. spring-wire guide

Percutaneous Sheath Introducer Procedural Kits With Silk Sutures

Contents:

- Arrow-Flex® Sheath: 9 Fr. (2.9 mm tip ID) x 10 cm and Tissue Dilator
 Spring-Wire Guide, Marked: .025" (0.64 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End "J" Tip on Other)
 with Arrow Advancer
 Catheter: 20 Ga. x 1-3/4" (4.45 cm) Radiopaque over 22 Ga. TW Introducer Needle
 Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) RW
 Pressure Transduction Probe
 Introducer Needle: 20 Ga. x 1-1/2" (3.81 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
 Injection Needle: SafetyGlide™1 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock Syringe
 Injection Needle: SafetyGlide™1 23 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
 Obturator: 8 Fr.
 Cath-Gard® Catheter Contamination Shield: 80 cm with TwistLock™ Adapter
 Dust Cap: Non-Vented
 SharpsAway® Il Locking Disposal Cup
 SharpsAway® Disposal Cup
 Maximal Barrier Drape™ with 4" fenestration
 Drape: 26" x 18" (66 cm x 46 cm) with adhesive
 Safety Scalpel: #11
 Stopcock: 4-way High-Flow
 Prep Sponge Swab
 Gauze Pad: 2" x 2" (5 cm x 5 cm)

9 Fr.

- Gauze Pad: 2" x 2" (5 cm x 5 cm)


 Gauze Pad: 4" x 4" (10 cm x 10 cm)

 Suture: 3-0 Braided Silk with Straight Needle

¹A trademark of Becton, Dickinson and Company

Rx only
Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.
California Prop. 65

MARNING: Cancer and Reproductive Harm.
www.P65Warnings.ca.gov
Not made with natural rubber latex.
Sterilized using ethylene oxide unless otherwise indicated in the contents list.
Fluid path components are non-pyrogenic.























(01)XXXXXXXXXXXXXXX (17)YYMMDD (11)YYMMDD

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex®

LBL075685 R00 (2024-05)















