



REF NA-04020-C (IPNXXXXXX)

20<sup>Ga.</sup> 4.45<sup>cm</sup> catheter length .018<sup>inch dia.</sup> spring-wire guide

# Arterial Catheterization

**Contents:**

- 1: Catheter: 20 Ga. x 4.45 cm Radiopaque Polyurethane over 22 Ga. TW Introducer Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide
  - 1: Injection Needle: SafetyGlide™<sup>1</sup> 25 Ga. x 1" (2.54 cm)
  - 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
  - 1: Syringe: 3 mL Luer-Lock
  - 1: 10 mL Luer-Lock Pre-Filled Saline Syringe
  - 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep
  - 2: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®<sup>2</sup> One-Step Solution with Hi-Lite Orange™ Tint
  - 1: Wing Clip: Suture
  - 1: SharpsAway® Disposal Cup
  - 1: SharpsAway® II Locking Disposal Cup
  - 1: Drape: 36" x 41" (91 cm x 104 cm) with 3" x 5" (8 cm x 13 cm) fenestration
  - 1: Needle Holder
  - 1: Safety Scalpel: #11
  - 1: Checklist Sheet
  - 1: Medication Label: 1% Lidocaine
  - 4: Gauze Pad: 2" x 2" (5 cm x 5 cm)
  - 4: Gauze Pad: 4" x 4" (10 cm x 10 cm)
  - 1: Dressing: CHG 7 cm x 8.5 cm
  - 1: Transducer Cover: 14 cm x 147 cm
  - 1: Surgical Apparel: Bouffant Cap
  - 1: Surgical Apparel: Mask with Eye Shield
  - 1: Suture: 3-0 Silk with Curved Needle
  - 1: Tubing±: 6" Extension with One-Way Stopcock
- <sup>1</sup>A trademark of Becton, Dickinson and Company.
- <sup>2</sup>A registered trademark of CareFusion or one of its subsidiaries.

**Rx only**

**Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. [www.teleflex.com/IFU](http://www.teleflex.com/IFU)**

**California Prop. 65**

**⚠ WARNING:** Cancer and Reproductive Harm. [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)

**±Precaution: Contains Phthalate: DEHP.**

**Warning: The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction. Proceeding from the present state of scientific knowledge, risks for male premature infants cannot be excluded in the case of long-term exposure or application. Medical products containing phthalates should be used only temporarily with pregnant women, nursing mothers, babies and infants.**

**Not made with natural rubber latex.**

Store between 20 – 25°C (68 – 77°F).

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.

Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.



**LOT** Lot Number

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXXXX  
(17)YYMMDD  
(11)YYMMDD  
(10)Lot Number

Product of xxx  
Packaged in xxx

Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA

**Teleflex®**  
[www.teleflex.com](http://www.teleflex.com)

LBL076358 R00 (2024-09)

Arterial REF NA-04020-C

20<sup>Ga.</sup> 4.45<sup>cm</sup> catheter length .018<sup>inch dia.</sup> spring-wire guide



# Arterial Catheterization

