



ARROW®

REF NA-04020-K (IPNXXXXXX)

20 Ga. 4.45 cm catheter length .018 inch dia. spring-wire guide

Radial Artery Catheterization


Contents:

- 1: Catheter: 20 Ga. x 4.45 cm Radiopaque Polyurethane over 22 Ga. TW Introducer Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX
- 1: Injection Needle: Eclipse™¹ 25 Ga. x 1" (2.54 cm) Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX
- 1: Filter Needle: 18 Ga. TW x 1-1/2" (3.81 cm) with Blunt Fill Tip Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX
- 1: Syringe: 3 mL Luer-Lock Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX
- 1: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®² One-Step Solution with Hi-Lite Orange™ Tint Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX
- 1: Wing Clip: Suture Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX
- 1: SharpsAway® II Locking Disposal Cup Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX
- 1: Towel: 13" x 18" (33 cm x 46 cm) Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX
- 1: Drape: 18" x 26" with 2" x 3" fenestration, with adhesive Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX
- 1: Scalpel: #11 Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX
- 4: Gauze Pad: 2" x 2" (5 cm x 5 cm) Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX
- ¹A trademark of Becton, Dickinson and Company.
- ²A registered trademark of CareFusion or one of its subsidiaries.

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU

California Prop. 65
⚠ **WARNING:** Cancer and Reproductive Harm.
www.P65Warnings.ca.gov

Not made with natural rubber latex.
Store between 15 – 30°C (59 – 86°F). Avoid freezing and excessive heat above 40°C (104°F). Sterilized using ethylene oxide unless otherwise indicated in the contents list. Fluid path components are non-pyrogenic.



LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Set Product of xxx
Packaged in xxx

Teleflex®
www.teleflex.com

LBL077607 R00 (2025-05)



ARROW®

REF NA-04020-K (IPNXXXXXX)

20 Ga. 4.45 cm catheter length .018 inch dia. spring-wire guide

Radial Artery Catheterization


Contents:

- 1: Catheter: 20 Ga. x 4.45 cm Radiopaque Polyurethane over 22 Ga. TW Introducer Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX
- 1: Injection Needle: Eclipse™¹ 25 Ga. x 1" (2.54 cm) Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX
- 1: Filter Needle: 18 Ga. TW x 1-1/2" (3.81 cm) with Blunt Fill Tip Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX
- 1: Syringe: 3 mL Luer-Lock Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX
- 1: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®² One-Step Solution with Hi-Lite Orange™ Tint Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX
- 1: Wing Clip: Suture Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX
- 1: SharpsAway® II Locking Disposal Cup Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX
- 1: Towel: 13" x 18" (33 cm x 46 cm) Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX
- 1: Drape: 18" x 26" with 2" x 3" fenestration, with adhesive Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX
- 1: Scalpel: #11 Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX
- 4: Gauze Pad: 2" x 2" (5 cm x 5 cm) Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX
- ¹A trademark of Becton, Dickinson and Company.
- ²A registered trademark of CareFusion or one of its subsidiaries.

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU

California Prop. 65
⚠ **WARNING:** Cancer and Reproductive Harm.
www.P65Warnings.ca.gov

Not made with natural rubber latex.
Store between 15 – 30°C (59 – 86°F). Avoid freezing and excessive heat above 40°C (104°F). Sterilized using ethylene oxide unless otherwise indicated in the contents list. Fluid path components are non-pyrogenic.



LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Set Product of xxx
Packaged in xxx

Teleflex®
www.teleflex.com

LBL077607 R00 (2025-05)



REF NA-04020-K

ARROW®

20 Ga. 4.45 cm catheter length .018 inch dia. spring-wire guide

Radial Artery Catheterization



REF NA-04020-K

ARROW®

20 Ga. 4.45 cm catheter length .018 inch dia. spring-wire guide

Radial Artery Catheterization