



REF CDC-15123-XPN1A

(IPNXXXXXX)

12^{Fr.} 3^{Lumen} 20^{cm catheter length} .035^{inch dia. spring-wire guide}

ARROW®

ErgoPack[®]
Complete

Pressure Injectable Three-Lumen Hemodialysis Catheterization for High Volume Infusions



Contents:

- | | |
|---|---|
| 1: Three-Lumen Catheter: 12 Fr. (4.06 mm OD) x 20 cm Pressure Injectable Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps | 1: Safety Scalpel: #11 |
| 1: Spring-Wire Guide, Nitinol, Marked: .035" (0.89 mm) dia. x 27-1/2" (70 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow® GlideWheel™ Wire Advancer | 1: Checklist/CLIP Sheet |
| 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle | 1: Sterile Procedure Sign |
| 1: Injection Needle: Eclipse™ ¹ 25 Ga. x 1" (2.54 cm) | 1: Medication Label: 1% Lidocaine |
| 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm) | 2: Gauze Pad: 2" x 2" (5 cm x 5 cm) |
| 1: Pressure Transduction Probe | 5: Gauze Pad: 4" x 4" (10 cm x 10 cm) |
| 1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe | 1: Surgical Apparel: Impervious Gown |
| 1: Injection Needle: SafetyGlide™ ¹ 23 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe | 1: Dressing: I.V. Advanced 8.5 cm x 11.5 cm |
| 1: Syringe: 3 mL Luer-Lock | 1: Dressing: BIOPATCH® ³ |
| 1: Tissue Dilator: 10 Fr. (3.5 mm) x 10.2 cm | 1: Transducer Cover: 14 cm x 147 cm |
| 1: Tissue Dilator: 12 Fr. (4.0 mm) x 14 cm | 1: Surgical Apparel: Mask with Eye Shield |
| 3: 10 mL Luer-Lock Pre-Filled Saline Syringe | 1: Surgical Apparel: Bouffant Cap |
| 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep | 1: Suture: 3-0 Silk with Curved Needle |
| 2: 3 mL Applicator 2% CHG and 70% IPA Chloraprep® ² One-Step Solution with Hi-Lite Orange™ Tint | 1: Tubing±: 20" (50 cm) Extension |
| 1: SharpsAway® II Locking Disposal Cup | 1: HemoHopper® Fluid Receptacle |
| 1: SharpsAway® Disposal Cup | |
| 1: Maximal Barrier Drape™ with 4" fenestration | |
| 1: Towel | |
| 1: Needle Holder | |

¹A trademark of Becton, Dickinson and Company.
²A registered trademark of CareFusion or one of its subsidiaries.
³A registered trademark of Johnson & Johnson Corporation.

±Precaution: Contains Phthalate: DEHP.

Warning: The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction. Proceeding from the present state of scientific knowledge, risks for male premature infants cannot be excluded in the case of long-term exposure or application. Medical products containing phthalates should be used only temporarily with pregnant women, nursing mothers, babies and infants.

Not made with natural rubber latex.

Store between 20 – 25°C (68 – 77°F).

Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.

Contraindications: The Multi-Lumen Acute Hemodialysis Catheter (AHDC) is not designed for long-term (≥ 30 days) hemodialysis or for insertion into thrombosed vessels.

